

Product datasheet for **SC209079**

SLC2A11 (NM_001024939) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	SLC2A11
Synonyms:	GLUT10; GLUT11
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001024939
Insert Size:	1493 bp



Insert Sequence:	>SC209079 3'UTR clone of NM_001024939 The sequence shown below is from the reference sequence of NM_001024939. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CTGGAGGTTATCCAGTCAACAGAACTCTAGTCCCAAAGGGGTGGCCAGAGCCAAAGCCAGCTACTGTCC TGTCCTCTGCTTCTGCGAGGCGCCCTGGTCTCACTCCCTCTGCATTCTCATTAAAGGAGTGTTAT TGAGCACCCCTTTGTGTGCAGACATGGCTCCAGGTGCTTAGCAATCAATGGTGAGCGTGGTATTCCAGGC TAAAGGTAATTAAGTACAGAAAAATCAGTAACAACATAATTACAGGCTGGTTGTGGCAGCTCATGACTG TAATCCCAGCACTTTGGGAGGCCAAGGTGGGAGGATCAATTGAGGCCAGAGTTGAAACCGAGCTAGGT AACATAGTGAGACCCCTATCTCTACAAAAAATTTAAACATTAGCTGGGCATGGTGGTATGTGCTAAC AGCTCTAGCTACTCAGGAGGCTGAGGCAGCAGGATCACTTGAGTCCAAGAGTTCAAGGTAGCAGTAAGC TACAATCACACCACTGCATGCCAGACTGGGTGACAGAGGGAGACTTCATCTCTTAAACATAATAATA ATAATTACAGACTCAGGAAATGCAGTGAAAGAAAAATACAGTTGGCCAGGTGAGGTGGCTGATGCCTG TAATCCCAGCACTTTGGGAGGCCAAGATGGGAAGATTGCTTTGAGACCAGAAGTTTGAGACCAGCCTGG GCCACATAGTAAGATCCTGTTTCTACCAAAAAAAAAAAAAAAAAAATTAGCTGGGTGTGGTGGTACAT GCCTGTGGTCCCAGCTACTCAGGAGGCTGAAATGGGAGGATCACTTGAGCCTGGGAGGTCGAGGCTGCA GTGAGTCCTAATTGAGCACTGCACTCCAGCCTGGACAACAGAAAGAGACCATGTTTCAAAAAAATAA ATACAGGTTGTAGTGGGTATGGGTATGCATACCAGGAGGCCGACCTCGTCTGAGAGGGGAGTGGTCAG AAAAGAATTTTCTGAGGAATTGATGTTTTTAATCAACTTTATTGAGGTATAATTTATACATAACCAACT GCATCCATTTTAAGTATATATTTGGCGAGTTTTGAGTTCTAAGTATAGTTTGGTGAGTGATACACTT GGGAAACACCACCGTGATCAAGATGGAACCTTTACCCTACCCCAAGGCACCCACATGCCCATTTGCTAT CAGTAACCCACCCAGTCCCAGACCTGGGGAACCACTGACCTGCTTTCTGCCTCAATTGGTCAGATTTG CCTTTTTTCAGAATCTGAAGTCATTCGTAAGTATGTACACATTTGTGTCTGGCTTGGCTCCAGATAA AGTTTTTGGGATTCATGAATGTTGTGCACGTATTAGGAGAGACTCCTTTGTATTGCTGAGTAGTATTC CCCTCTGTGGTTAGACCATGATTTATTTATCCATCTACCTGTTGGTGAACATTTGGCTGTTTCTAATT CTTGCCATCATGAATAAACTGCTGTGAATGTTCTACAATAA ACGCGTAAGCGGCCCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001024939.4
Summary:	This gene belongs to a family of proteins that mediate the transport of sugars across the cell membrane. The encoded protein transports glucose and fructose. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Locus ID:	66035

MW: 55.7