

Product datasheet for **SC209081**

UBE2L6 (NM_004223) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: UBE2L6 (NM_004223) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: UBE2L6
Synonyms: RIG-B; UBCH8
ACCN: NM_004223
Insert Size: 726 bp
Insert Sequence: >SC209081 3'UTR clone of NM_004223
The sequence shown below is from the reference sequence of NM_004223. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCCGATTGGAGTGGACCGCCCTCCTAACTCATGTTCTGACCCTCTGTGCACTGGATCCCTCGGCATA
GCGGACGGACACACCTCATGGACTGAGGCCAGAGCCCCCTGTGGCCATTCCCCATTCATTTTTCCCTT
CTTAGGTTGTTAGTCATTAGTTTGTGTGTGTGTGGTGGAGGGAAGGGAGCTATGAGTGTGTGTGTG
TGTATGGACTCACTCCCAGGTTACCTGGCCACAGGTGCACCCTTCCCACACCCTTTACATTCCCCAGA
GCCAAGGGAGTTTAAGTTTGCAGTTACAGGCCAGTTCTCCAGCTCTCCATCTTAGAGAGACAGGTCACC
TTGCAGGCCTGCTTGCAGGAAATGAATCCAGCAGCCAACCGAATCCCCCTAGGGCTCAGGCACTGAGG
GCCTGGGACAGTGGAGCATATGGGTGGGAGACAGATGGAGGGTACCCTATTTACAACAGTGCAGCCA
AGCCACTGATGGGAATATACAGATTTAGGTGCTAAACATTTATTTCCACGGATGAGTCACAATCTGA
AGAATCAAACCTCCATCCTGAAAATCTATATGTTTCAAACCACTTGCCATCCTGTTAGATTGCCAGTT
CCTGGGACCAGGCTCAGACTGTGAAGTATATATCCTCCAGCATTAGTCCAGGGGGAGCCACGGAAC
CATGTTCTTGCTTAAGCCATTAAGTCAGAGATGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_004223.5
Summary:	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is highly similar in primary structure to the enzyme encoded by the UBE2L3 gene. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2011]
Locus ID:	9246
MW:	27.1