

Product datasheet for **SC209082**

ASAHL (NAAA) (NM_001042402) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ASAHL

Synonyms: ASAHL; PLT

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_001042402

Insert Size: 726 bp

Insert Sequence: >SC209082 3'UTR clone of NM_001042402
The sequence shown below is from the reference sequence of NM_001042402. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AACCTCAGCCTGGAGGCACTTTTCCAGTGAAGAAGAACCATCAGGTGATTGTCAGGCCTCTGAGC
CCAAGCTAAGCCATCATATGCCCTGTGACCTGCACGTATACATCCAGATGGCCTGAAGCAACTGAAGAT
CCACAAAAGAAGAGAAAATAGCCAGTTCCTGCCTTAAGTATGACATTTACCATTTGTGATTTGTTCTT
GCCCCACCCTAACTGATCAATTGACTTTGTGACAAGACACCCTCCCCACCCTTGCAATAATGTACTTTG
TGATATTTCCCCGCCCTTGTAATGTACTTTGTACGATACACCCTCCCCACCCTTGAGAAGGTACTTTG
TAATATCCTCTCGCCCTTGAGAATGTACTTTGTAAGATCCACTTCCTGCCTGCAAAAAATTGCTCCGA
ACTCCACCGTCTATCCCAAACCTATAAGAACTAATGATAATCCTACCACCCTTTGCTGACTCCTTTTTT
GGACTCAGCCACGTGCACCCAGGTGATTAAGGCTTTATTGCTCAGATGAAGCCTATTTGGTGGTCT
CTTCACACAGATGTGTGTGACAGTGATAACAGTGTCACGTGGTTTAACTACAATGTGTCCAGATATTGA
ATTTTGTGTTAAGCCTGCTTGGGATTTGTTGTGCTCCCAAAATTTCAAGTTTGTGTCTTTTGTCAATTT
TGGAAAATAGTCAGCTATCTCAGAATTGCTCACCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001042402.2</u>
Summary:	This gene encodes an N-acyl ethanolamine-hydrolyzing enzyme which is highly similar to acid ceramidase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	27163
MW:	27.6