

Product datasheet for **SC209109**

MSANTD2 (NM_024631) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: MSANTD2

Synonyms: C11orf61

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_024631

Insert Size: 708 bp

Insert Sequence: >SC209109 3'UTR clone of NM_024631
The sequence shown below is from the reference sequence of NM_024631. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATTGGATACCCCTTAAATTTAACAACATGATGTCATCCTTCATAAGGATTTGGGCTCTTAGCTCCTTC
TTCTCTACTCACTTCCCATACCCGGACCACCCCTCATCCAGATGCCGCCGTCAGACTCTTCATGGAAA
CCCTTTCTTCAATTGGGCCAGTACGACTTCTAAGGAATCATAGTAGACTTGGGCCAGAAAAACAGACC
TCACCTGAAAATGCTCATTTTGGCAAGCAAGATATATAGCGAACTTGCATGGTGGGTCAGCTATTTCT
TTTTTAAAAACAGAAAAACAATTTTAAATGGATTTTAGAGCCCTAATATAAACAAATGTAGGTACATT
CCATATTGGACAAAACCTATACTTTGAGTTTCATATGAAATTGAAAAATTGTATTTTTTTCACAATGTT
TCATTTTACACATCTCTATGTCTCTATATAAGTAATATTTACAACCTGAATAAATAAGCAATAGATA
ATGGCAAGTTAAATCTTGATTACAGTTAAAAAGACTTTTCAATATCACCCCTTAGCTACTATTGTATA
TAATTTAGATTTTCATTTTTCACCAGCCAAGTGTGTTGCTATTTCCAATCAGTATTGCTTGCAAGA
CCTTTGTTTCTGTGATGTACCAACTTTACAGTTGTGTGCCATTTAAACCTTTTTTAAAAA
TTAAAGTAATTCCTGGCC
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024631.4</u>
Locus ID:	79684
MW:	27.6