

Product datasheet for SC209117

prothymosin, alpha (PTMA) (NM_001099285) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	prothymosin, alpha (PTMA) (NM_001099285) Human 3' UTR Clone
Symbol:	prothymosin, alpha
Synonyms:	TMSA
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001099285
Insert Size:	730 bp
Insert Sequence:	>SC209117 3'UTR clone of NM_001099285 The sequence shown below is from the reference sequence of NM_001099285. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGAAGCAGAAGACCGACGAGGATGACTAGACAGCAAAAAAGGAAAAGTTAACTAAAAAAAAAAGGC
CGCCGTGACCTATTCACCTCCACTTCCCGTCTCAGAATCTAAACGTGGTCACCTTCGAGTAGAGAGGC
CGCCCCGCCACCGTGGCAGTGCCACCCGAGATGACACGCGCTCTCCACCACCAACCCAAACCATG
AGAATTTGCAACAGGGGAGGAAAAAGAACCAAACTTCCAAGGCCCTGCTTTTTTCTTAAAAGTACT
TTAAAAAGGAAATTTGTTTGTATTTTTTTTACATTTTATATTTTTGTACATATTGTTAGGGTCAGCC
ATTTTTAATGATCTCGGATGACCAAAACAGCCTTCGGAGCGTTCTCTGTCTACTTCTGACTTTACTTG
TGGTGTGACCATGTTTATTATAATCTCAAAGGAGAAAAAACCTTGTAAGCAAAATGACAAC
AGAAAAACAATCTTATCCGAGCATTCCAGTAACCTTTTTTGTGTATGTACTTAGCTGTACTATAAGTAG
TTGTTTGTATGAGATGGTTAAAAAGGCCAAAGATAAAAGGTTTCTTTTTTTTCCTTTTTGTCTATG
AAGTTGCTGTTTATTTTTTTGGCCTGTTTGTATGTGTGAAACAATGTTGTCCAACAATAAACAGG
AATTTTATTTTGTGAGTTGTTCTAACAAAGCTGTCTCAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001099285.2
Summary:	Prothymosin alpha may mediate immune function by conferring resistance to certain opportunistic infections.[UniProtKB/Swiss-Prot Function]
Locus ID:	5757
MW:	28.7