

Product datasheet for **SC209151**

COMMD7 (NM_053041) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	COMMD7 (NM_053041) Human 3' UTR Clone
Symbol:	COMMD7
Synonyms:	C20orf92; dj1085F17.3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_053041
Insert Size:	713 bp
Insert Sequence:	>SC209151 3'UTR clone of NM_053041 The sequence shown below is from the reference sequence of NM_053041. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTCAGAACCATGAGTGTCTCTGCTATTTCTGTCCCTGCATCTCCCCTGGCCCCGTTCCCTGCCC
TCCTCCCTTCCCTGGGTGACTGCTCTGAGAGGCACTTCACTCACAGGCCTGTGGGATGCTCCATGGGGC
CCTGCTGGCTCCATGGGCCCCAGGTGCAAAGGGTTTCTGAAAAACAGCAGGATTAAGTACTGAAAGAGC
CCAACACAATTACCCTGTAACTCTCTGTTAGGGCAACCACCACCTGTCTTCCAGGACACATTTTT
AGATACTCTGACAGGCCACTGCATCTCAGATTCAGGGGAGAAAAATAAGTTGTACCTCCCTTCAAAGT
TCCAGAGTAAACAAATGGTGCCATCATTCAAGATAACATGCTGATCACCTCCTCCCAAAAAGCAAGAG
CTTGTTTATGGCTGAGGAATCGGCGGATTGTCTGAATGACACATATACAGAGCCCCACGGATTTCTGC
ACACTCTGGGTCTGTGCTGGTGGAACTTGCCAATCAGTTCTTAATGAGGCACCTGTGTGTAATACAT
GCTTGGTCTTCTCTGCAGAGAACTGAGGCTAACTCTGTCCCTACTTCTGGTTTTGCCCTGTCATGTCG
TAACGAGGTGGGCCTTTTGAGGCCATTTAGTTTGAGTTTCAGCCAACCACCTCTGTTGGTTAGATGAT
GAATAAAAAGTTCTGAAGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_053041.3</u>
Summary:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). Associates with the NF-kappa-B complex and suppresses its transcriptional activity (PubMed:15799966).[UniProtKB/Swiss-Prot Function]
Locus ID:	149951
MW:	26.7