

Product datasheet for SC209162

POTEF (NM_001099771) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	POTEF
Synonyms:	A26C1B; ACTB; POTE2alpha; POTEACTIN
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001099771
Insert Size:	739 bp
Insert Sequence:	>SC209162 3'UTR clone of NM_001099771 The sequence shown below is from the reference sequence of NM_001099771. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCCTCCATTGTCCACCGCAAAATGCTTGAGGTGGACTCTGACTTAGTTGCGTTACACCCTTTCTTGACA AAACCAAACTTCTCAGAAAACAACATGAGATTGGCATGGCTTTATTTGTTTCTTGTTTCATTTTTTGT TTTGTTTTTATTGGCTTGACTCAGGATTTAAAAACCGGAATGGTGAAGGTGACAGCAGTCGGTTGGAG GAAGCTTCCTCCAAAGTTCTACAATGTGGCCAAGGACTTTGATTGTACCTTGTTCTTTTCAATAGT CATTCCAAATATTGTGAGACGATTGTTTCAGGAAGCCCTTGCCCTGCTAAAAGCCACCCCACTTCTC TCTAAGGAGAATGGCCAGTCCTCTCCCTAGTTCACACAGGGGAGGTGATAGCATTGCTTTTGTGCAAA TTACATAATGCAAAATTTTTGAATCTTCGCCTTAATACTTTTAATTTTGTGTTTATTTGAATGATCA GCCTTCGTGGCCCCCTCTTTTGTACCCCAACTTGGGGTGTATGAAGGCTTTTGGTCTCCCTGAGAGTG GCTGGAGGCAGCCAGGGCTTACCTGTACTCTGACTTGAGGAGAGTTGGATAAAAGTGACACCTTAAAA CAAATTGAGGAAGCACAGTATTTAGTACAGTGACAGCTTAGCATGTTGACAAGTGAATAAAATGC TCAGTTCTGAACTGGACAATGTAAGACACAACGAGGAAACTGGAAT ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001099771.2</u>
Locus ID:	728378
MW:	28.1