

Product datasheet for **SC209207**

ADAMDEC1 (NM_001145271) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ADAMDEC1 (NM_001145271) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ADAMDEC1
Synonyms: M12.219
ACCN: NM_001145271
Insert Size: 745 bp
Insert Sequence: >SC209207 3'UTR clone of NM_001145271
The sequence shown below is from the reference sequence of NM_001145271. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAGATGCTCAAACCATACCACAGAGTGAATCCAAAAGTCTGCTTCACTGAGATGCTACCTTGCCAGG
ACAAGAACCAAGAACTCTAACTGTCCAGGAATCTTGTGAATTTTCAACCATAATGGTCTTTCACCTGT
CATTCTACTTTCTATATTGTTATCAGTCCAGGAAACAGGTAAACAGATGTAATTAGAGACATTGGCTCT
TTGTTTAGGCCTAATCTTTCTTTTTACTTTTTTTTTCTTTTTCTTTTTTTAAAGATCATGAATTT
GTGACTTAGTTCTGCCCTTTGGAGAACAAAAGAAAGCAGTCTTCCATCAAATCACCTTAAATGCACGG
CTAAACTATTCAGAGTTAACTCCAGAATTGTTAAATTACAAGTACTATGCTTTAATGCTTCTTTCAT
CTTACTAGTATGGCCTATAAAAAAATAATACCACTTGATGGGTGAAGGCTTTGGCAATAGAAAGAAGA
ATAGAATTCAGGTTTTATGTTATTCCTCTGTGTTCACTTCGCCTTGCTCTTGAAAGTGCAGTATTTTC
TACATCATGTCAAGAATGATTCAATGTAATATTTTTTCATTTTATCATGTATATCCTATACACATCT
CCTTCATCATCATATATGAAGTTTATTTTGGAGAAGTCTACATTGCTTACATTTTAATTGAGCCAGCAA
GAAGGCTTAATGATTTATTGAACCATAATGTCAATAAAAAACAACTTTTGAGGC
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001145271.2</u>
Summary:	This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells. [provided by RefSeq, Jul 2008]
Locus ID:	27299
MW:	29.3