

## Product datasheet for SC209212

### ABHD3 (NM\_138340) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	ABHD3 (NM_138340) Human 3' UTR Clone
Symbol:	ABHD3
Synonyms:	LABH3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138340
Insert Size:	720 bp
Insert Sequence:	>SC209212 3'UTR clone of NM_138340 The sequence shown below is from the reference sequence of NM_138340. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATGTTGAGCATGGACATGAACCTCTCTTAACCATGTAGTCTTTGGGTGCATTTTGTCTGAACCACAATT
GTGAAGGCAGCTCAGCTTAGTGACAAAATTTAACTGTTGTATATAAAGCAAATAAGCCAGCAGATGGG
TGAAGAGGTCCAGAATGATATGCAAAAACACTTTTTAGAGAAAACAAAACAACCTTTGTAGCAACAAATT
AAATATAGTATTAGATTGTTACTTACGTAGATTTTATTTTACTATGCCTTACCAAGTACATCCTTAAA
CAAAGTAGTATGTACATGAAATTGCACTTAACCAAAACTATTGTGTAAAACAAATTTAATTCCTCAGG
GTTTTAATTTAACTAGTATTTTTTTAGATTATTTGTTTTAGGTGATTTAATGGTACTTTAATAACTAC
TAAGAAATATTGGCTATTTCAATGTAAGTTATAAGGTGGTACATTCCTAAGGGTATTTATAGTTGATGA
TAACATGAAAACGAAATAAGATAAAAACAACGTGCTAAATCTTTTATGTATTCTAACTTTAAAAGAC
AAGTGCAACAAAGTTAGACTGACTTCTATATGTGCTCTTTACTCTGATAATATTAATAGGACTAAC
TTATGTTTTATAATGATTATAATTTACATGCTTATTTTTAAAATAGTATATGTGGACACATATATATCA
TTATATTTAAAATAAATTCTACCATTTTTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_138340.5</a>
<b>Summary:</b>	This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined. [provided by RefSeq, Jul 2008]
<b>Locus ID:</b>	171586
<b>MW:</b>	28.5