

Product datasheet for SC209225

LYRM5 (ETFRF1) (NM_001001660) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LYRM5 (ETFRF1) (NM_001001660) Human 3' UTR Clone
Symbol:	LYRM5
Synonyms:	LYRM5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001660
Insert Size:	861 bp
Insert Sequence:	>SC209225 3'UTR clone of NM_001001660 The sequence shown below is from the reference sequence of NM_001001660. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TATTATTCAGATACCAACAAAATAATATCATTACTACTTTAATTTAGCTATCAGTGCCAGCTGTTT
ATGTATACCAGATGTTGTAATAAATTCTAACTTAAATGGGAAGATATACATGTTGTGTAATAAATCC
CTGAGCTGCCCTACTGAATAAATAGGTTTCACTTCTGTTTATACGGAGAAAGTATCAGCAACTTTAT
GCTCAATTTTGATACAAACATAGCAATTTAGCTATACTACCGATCATAAATTAATGAGCACCAATTTT
GAATGAAATATAATACTTAATCCTCTAACTTAATAGGACTTAGCCAATTATTTTAGCTTATTTCT
CCATTTTACTGACAAGAGAATGCCATTTACCGAACATTTAATAAGCAAAATAATTTCTTACACACTTAC
TTCTCTTAACTTTTGAAAGGAATCCCTTATTTTTTACCCAATTTGGTATATTAACAGCACATACTCAA
GGTTCACTACAAAACAAACAGTTCCTGGTAATGATTTAAATGTAGTTATAGAAATAAATAATATGTATG
GAGTCATTACTTCTGACCTTGAAATAGCCTGCTGGTGACTGGCATTAACTACATAACTAACTATTCAA
TTATTTCCATTATTAATGTTGCTGCTGCTACCTTTGGGCTAAACCAAGTCACCTGTTGTGTATCAAT
TACCTTCTTTGATAAAGGAAATAAATAAATGTTATCCTATTTTTTGTGATAAAGGCAGATTTGT
TTTGCATCTACTTCTAAATATGGTTTTAAGTTTAAGCAACACGCCTTTACATAATATCCACAGCACCT
TTTAGAAATAAAGGATACTTCTAACAAGCTGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001001660.3
Summary:	Acts as a regulator of the electron transfer flavoprotein by promoting the removal of flavin from the ETF holoenzyme (composed of ETFA and ETFB).[UniProtKB/Swiss-Prot Function]
Locus ID:	144363
MW:	33.8