

## Product datasheet for **SC209319**

### ZNF114 (NM\_153608) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | ZNF114 (NM_153608) Human 3' UTR Clone                                                                                                                                                                                                              |
| Symbol:                   | ZNF114                                                                                                                                                                                                                                             |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_153608                                                                                                                                                                                                                                          |
| Insert Size:              | 760 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC209319 3'UTR clone of NM_153608</p> <p>The sequence shown below is from the reference sequence of NM_153608. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACTCACAAGATGAGAAGCCCTGTGAATGAAGGAAGGTGGAAAAATTTTTCATTAATTTTCTGACTGTA
CCAAACATGTGAGGAGGACATATTGGAAGGGAGCTCAAGGGGTAGCATGAGTGAGAACATCTTCCCTG
AACTCTCGTATCTTACAGAAATGTGAAAAAAACCCTGTGAAGGTAAAGTCTACAGAAAGCCTTTCATC
TTCATTCATCTTGAGTAGACATTTGTTCTCACCTGGAGAGAACTGCGAATCTAAAAGGAATATGACA
AAGCCTTCAGCGTGGTCTCAAATTCATGGTTCATACAAGAACTCAGCTGCAGAGACTCCTTACGGAAA
TAAAAAATGTAGGAAAGACCTGCCGGCCGCGGTGGCTCATGCCTGTAATCCCAGCACTTTGGGAGGCCG
AGGCGGGCGGATCACGAGGTCAGGAGATCAAGACCATCCTGGCTAACACGGTGATACCCCGTCTCTACT
AAAAATACAAAAAAATTAGCCAGGCGTGGTGGCAGGTGCCTGTAGTCCCAGCTACTTGGGAGGCTGAG
GTAGGAGAATGGTGTAACCTGGGAGGTGGAGCTTGCACTGAGCCGAGATGGAGCCACTGCACTCCAGC
CTGGGGGAAAGAGCGAGACTCCGTCTCAAAAAAAAAAAAAAGTAGGAAAGACCTCTTTAAACACTTACC
ACTCATGTTGCATGTGAATCCACATCCTGAAGGGAGGCTACAAGTGAATAAAAATAAGAAAGCTTTCAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

|                    |                                                                                                                                                                                                                                                                                            |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | SgfI-MluI                                                                                                                                                                                                                                                                                  |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |


[View online »](#)

|             |                                                                                                                                                                                                |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq:     | <u>NM_153608.4</u>                                                                                                                                                                             |
| Summary:    | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]                                                                                                                  |
| Locus ID:   | 163071                                                                                                                                                                                         |
| MW:         | 29.2                                                                                                                                                                                           |