

Product datasheet for **SC209325**

AGTRAP (NM_001040197) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Product Name: AGTRAP (NM_001040197) Human 3' UTR Clone

Vector: pMirTarget (PS100062)

Symbol: AGTRAP

Synonyms: ATRAP

ACCN: NM_001040197

Insert Size: 740 bp

Insert Sequence: >SC209325 3'UTR clone of NM_001040197

The sequence shown below is from the reference sequence of NM_001040197. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTACCACATGTACCGGGAGCGCGGGGTGAGCTCCTGGTCCACACTGGTTTCCTTGGGTCTTCTCAGGA
CCGTAGTGCCTACCAGACGATTGACTCAGCAGAGGCCGCCAGATCCCTTTGCAGTCCCAGAGGGCAG
GAGTCAAGATGCCCGAGGGTACTGAAGCCAGCCAGCTGCGCCCCGCCCTGCCCGGGCCTTCTCGTG
CCTGGGAGGTCGTTCTAGGGATGCTCCTGACCTCCGTCTCTTGGACCTAAGATGGAATGTGTCCCCAGC
TCAGGGATTGCCTGAACCAAGAGGCCAGGAGCCCCATGGGCCGCCAGTACCATGCACACTCCTGTCC
CGAACTCCCTGAGGCCTCCCCTCCCTTCAGGGCACCCACTGGTTCCAGGCTGGAACCAGGGTCTCTCT
TTACCTCCTACCCCATGGTGGCACCACAGAGGCCCTCAGCCGAGTCTGCCTGAGTGTGCAAGCTCAG
GCCTTTAAGGACTGCTGATGCCCCCTCAGGCTCCCCAAGTTTGTGGGCTTTGGTGAAGCCCTGAG
AGCTTCAGGTCTGCTCAGCCCAGGAGCAGTCTGGCATGGGAGTGAAGCCCGTCTTCTCACTGCCT
GGTCACATGGTGCCTAGGGATGCAGGGCTGGAGGCCAGAGGTGTCAGCAACTGTGTCCCACCACAAC
CTCCAGCCTCCCTTTTTCAGAGCACAGCATTAAAGTTTGGGGAATTCTGTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001040197.2
Summary:	This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Locus ID:	57085
MW:	26.4