

Product datasheet for **SC209361**

ENSA (NM_207045) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ENSA (NM_207045) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ENSA
Synonyms:	ARPP-19e
ACCN:	NM_207045
Insert Size:	1628 bp



[View online »](#)

Insert Sequence:	<p>>SC209361 3'UTR clone of NM_207045 The sequence shown below is from the reference sequence of NM_207045. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AGCAAGCTTGCGGGTGGCCAAGTTGAAATGATGCTGCCCGGGCTCTGCCAGATCCTGAGACGCTTCCCC TCCTGCCACCCGGGTCTGTGCTGGCTCCTGCCCTTCTGTCTTTTGCAGCCAGGGGTGAGGAGGT GGCTCGGGTGTGGGCTGGAGAGGCAGAAGCCCTTCTGTTGGTGTCCAGCACATGGAGCCCTTGGG CTGAGCACCAAGACCTTGAACCTTTTTTGTTCCTTTTTTCAAATAACAGTTGGGAGAAATATCAA TGAAATTCTGGGGTGGGGTGGGGTGAAGGGTGGGGTGGGAGATATGGAGGAGTATGAATTAGGGC TTGGAGTTGGTAAAAACATTCTGACTATCCTTCTTAACACGTGGCTGATGTGGGTAGGTATGAGGG GAAGGAAGTGGAGTAGCCTAATGAAAAGGGTCTAGTTGAGCTCTGTAGATAAATGCCTTGTTCAGT GTGGTTGGAGACCTGGTGTGAGATAAAAGAACTCCATCCGCACAGACAGATGCAACAGCTCCTCTAG TTCTGCAGAGCTAGTTGAGAACTCAACATTAATCATTTTTAAAAAGTACTGTCCTTGAATAGATTGCT GTGGGAAGAAGGGCAGTGAGTGTGGGAGAAAGGAGCCGTGAGCGTGGGGAACCCACAGAGCCAAAGG ACTTTTTCAGTATTCGAAATAAACAAAAACAAACCCATGAAAAACCCACCCAGAAATATTCTGAAG CAGTTGGTGTGTTTTATTGGGGTTGGGAGGGGAAGAAAACAATACTGAAGACATTGTTAATGTGAGA TGCGGCTGAGTCCCAATACCATCACTTAGGAAAGGGAGCAGGAATGTTGGCATAAACCTTGCCATC CCTTTGTTCTGTAAACCAATGTAGAAGTGCATCAACCTCAGTTTGTGGCTCACAAAATGAGGCATTTA AGATATTTGATCTAAATGGTTTCTTGACACCTAAGCCTAGTCCAGTTTAGCAGAAAAGGAGGCCAGT CTCCTTTCCCTCTGCTTCCATTCCTTTCAAGTTCCCATGTTTGTGAGGGAGAAATGTGACAAACCA AAAGAAATATATTTCTTTGTAATAAGCCAAAAGAACTCAAGTCTTTCTGCTGGTGGCCAGCAGGA GGCAGAGTGATTTCTTCTCATGATGTGCCAATCTGCATTGCTGGAACCTGCTACTGAAAGAAGTTA AAAGCTGCCTGTGAATGTAACTGACAAAAGAAAGTTACCAAAGGAGAAAAAACTCCAGCGGGTGTACG GTGGCTCACGCCTATAATTCCAGCACTTTGGGAGGCCAAGGTGGGTGGATCACTTGAGGTGAGGAGTTT GAGACCAGCCTAGCCAACATGGTGAACCCCATCTCTACTAAAAATACAAAAATTAGCCGGGCGTGGCG GCATATGCCTATAATCCCACTTACTTGGGAGGCTGAGGTAGGAGAATCGTTGAACCCAGGAGGCAGAG TTTGAGTGAGCCGAGATTGTGCCATTGCACTCCAGCCTGGGCGACAGAGCGAGACTCCATCTAAAAAA AATAAATGAATAAAATAAAACAAGTACAATTTAACTAGA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_207045.2

Summary:

The protein encoded by this gene belongs to a highly conserved cAMP-regulated phosphoprotein (ARPP) family. This protein was identified as an endogenous ligand for the sulfonylurea receptor, ABCC8/SUR1. ABCC8 is the regulatory subunit of the ATP-sensitive potassium (KATP) channel, which is located on the plasma membrane of pancreatic beta cells and plays a key role in the control of insulin release from pancreatic beta cells. This protein is thought to be an endogenous regulator of KATP channels. In vitro studies have demonstrated that this protein modulates insulin secretion through the interaction with KATP channel, and this gene has been proposed as a candidate gene for type 2 diabetes. At least eight alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]

Locus ID:

2029

MW:

60.8