

Product datasheet for **SC209380**

DBNDD2 (NM_001048222) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: DBNDD2 (NM_001048222) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: DBNDD2
Synonyms: C20orf35; CK1BP; HSMNP1
ACCN: NM_001048222
Insert Size: 742 bp
Insert Sequence: >SC209380 3'UTR clone of NM_001048222
The sequence shown below is from the reference sequence of NM_001048222. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGTCAGCCAGAGCCGGATGTCAGGCTCTGAACGAGGCTACAAGGCTGGGCTGGGGAAGTACACAAGGA
TGGACAACCAATTTGGAGGAGCTGAGCCTGCCGGTGCTACATCAGACAGGACCACATCTAGGACCTCCT
CCTCCTCCTCCTCCGACTCCTCCACCAACCTGCATAGCCCAAATCCAAGTGATGATGGAGCAGATACGC
CCTTGGCACAGTCGGATGAAGAGGAGGAAAGGGGTGATGGAGGGGAGAGCCTGGAGCCTGAGCTAGC
AGTGGGCCCTGCCTACAGACTGACCACGCTGGCTATTCTCCACATGAGACCACAGGCCAGCCAGAGC
CTGTGGGAGAGAAGACCAGACTCTTTACTTGCAGTAGGCACCAGAGGTGGGAAGGATGGTGGGATTGTGT
ACCTTTCTAAGAATTAACCCTCTCCTGCTTTACTGCTAATTTTTTCTGCTGCAACCCTCCCACCAGTT
TTTGGCTTACTCCTGAGATATGATTTGCAAATGAGGAGAGAGAAGATGAGGTTGGACAAGATGCCACTG
CTTTTCTTAGCACTCTTCCCTCCCCTAAACCATCCCGTAGTCTTCTAATACAGTCTCTCAGACAAGTGT
CTCTAGATGGATGTGAACCTTAACTCATCAAGTAAGGTGGTACTCAAGCCATGCTGCCTCCTTACAT
CCTTTTTGGAACAGAGCACGGTATAAATAATAAACTAATAATAATATGCCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001048222.3
Summary:	May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro). [UniProtKB/Swiss-Prot Function]
Locus ID:	55861
MW:	28.3