

Product datasheet for SC209396

MCEMP1 (NM_174918) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MCEMP1 (NM_174918) Human 3' UTR Clone
Symbol:	MCEMP1
Synonyms:	C19orf59
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_174918
Insert Size:	745 bp
Insert Sequence:	>SC209396 3'UTR clone of NM_174918 The sequence shown below is from the reference sequence of NM_174918. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CAGAAAATGCCACAGTCTCACCTCAAATGAGAGGACATTGTGGCAGCCAAAGCCACAACCTGGAA
GATGGGGCTGCACCTGCCAACGAAGACGGGAAATGACCCCCCCCCCAGCCTAGTGTGAACCTGCCCC
TCGTCCCACGTATAGAAAAACCTCGAGTCATGGTGAATGAGTGTCTCGGAGTTGCTCGTGTGTGTAC
ACCTGCGTGCCTGTGTGTGCGTGTGTGCGCGTGTTCGTGTATGTGCGTGTGTGCGTGCCTGTGTG
TGCAATTTGCAAAGGGTGGACATTTCAAGTGTATCTCCAGAAAGGTGATGAATGAATAGGACTGAGAGT
CACAGTGAATGTGGCATGCATGCCTGTGTGTGTGACATATGTGAGTCTCGGCATGTACGGTGGGTGG
CTGTGTCTGAGCACCTCCAGCAGATGTCACTCTGAGTGTGGGTGTTGGTGACATGCATTGCACGGCCT
GTCTCCCTGTTTGTGTAACATACTAGAGTATACTGCGGCGTGTCTGTCTACCCATGTATGGTGG
GGGAGATTATCTCCGTACATGTGGGTGTGCCATGTGTGCCCTGTCACTATCTGTGGCTGGGTGAACG
GCTGTGTCAATTATGAGTGTGCCGAGTTATGCCACCCTGTGTGCTCAGGGCACATGCACAGACATTTA
TCTCTGCACTCACATTTTGTGACTTATGAAGATAAATAAAGTCAAGGGAAAAACAG
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_174918.3</u>
Summary:	This gene encodes a single-pass transmembrane protein. Based on its expression pattern, it is speculated to be involved in regulating mast cell differentiation or immune responses. [provided by RefSeq, Jul 2008]
Locus ID:	199675
MW:	27.3