

Product datasheet for **SC209429**

UGT1A7 (NM_019077) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: UGT1A7 (NM_019077) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: UGT1A7
Synonyms: GNT1; hUG-BR1; UDPGT; UDPGT 1-7; UGT-1A; UGT-1G; UGT1; UGT1-01; UGT1-07; UGT1.1; UGT1.7; UGT1A; UGT1A1; UGT1G
ACCN: NM_019077
Insert Size: 771 bp
Insert Sequence: >SC209429 3'UTR clone of NM_019077

The sequence shown below is from the reference sequence of NM_019077. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGAAAGCCACAAAATCCAAGACCCATTGAGAAAGTGGGTGGGAAATAAGGTAAAATTTGAACCATTCC
CTAGTCATTTCCAAACTTGAAAACAGAATCAGTGTAAATTCATTTTATTCTTATTAAGGAAATACTTT
GCATAAATTAATCAGCCCAGAGTGCTTTAAAAAATTCTCTTAAATAAAAATAATAGACTCGCTAGTCA
GTAAAGATATTTGAATATGTATCGTGCCCTCTGGTGTCTTTGATCAGGATGACATGTCCATTTTTC
AGAGGACGTGCAGACAGGCTGGCATTCTAGATTACTTTTCTTACTCTGAAACATGGCCTGTTGGGAGT
GCGGGATTCAAAGGTGGTCCCACGGCTGCCCTACTGCAATGGCAGTTTTAATCTTATCTTTGGCTT
CTGCAGATGGTTGCAATTGATCCTTAACCAATAATGGTCAGTCCATCTCTGTCGTGCTTCATAGGTG
CCACCTTGTGTGTTAAAGAAGGGAAGCTTTGTACCTTTAGAGGTAGGTGAAATGAATGAATGGCTTG
GAGTGCCTGAGAACAGCATATGATTTCTGCTTTGGGAAAAAGAATGATGCTATGAAATGGTGGGT
GGTGTATTTGAGAAGATAATCATTGCTTATGTCAAATGGAGCTGAATTTGATAAAAACCCAAAATACAG
CTATGAAGTGCTGGGCAAGTTACTTTTTTCTGATGTTTCTACAATAAAAATAAATAAATTT
ATATAAATTCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_019077.3</u>
Summary:	This gene encodes a UDP-glucuronosyltransferase, an enzyme of the glucuronidation pathway that transforms small lipophilic molecules, such as steroids, bilirubin, hormones, and drugs, into water-soluble, excretable metabolites. This gene is part of a complex locus that encodes several UDP-glucuronosyltransferases. The locus includes thirteen unique alternate first exons followed by four common exons. Four of the alternate first exons are considered pseudogenes. Each of the remaining nine 5' exons may be spliced to the four common exons, resulting in nine proteins with different N-termini and identical C-termini. Each first exon encodes the substrate binding site, and is regulated by its own promoter. The enzyme encoded by this gene has moderate glucuronidase activity with phenols. [provided by RefSeq, Jul 2008]
Locus ID:	54577
MW:	28.9