

Product datasheet for SC209431

KBTBD4 (NM_018095) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KBTBD4 (NM_018095) Human 3' UTR Clone
Symbol:	KBTBD4
Synonyms:	BKLHD4; HSPC252
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018095
Insert Size:	763 bp
Insert Sequence:	>SC209431 3'UTR clone of NM_018095 The sequence shown below is from the reference sequence of NM_018095. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTTACCAATTTGCAGTTTGTGTTGGCCTAAGGCTGTGGGAGGGGAGGAGAACTGCTCACTCCTTTTCC
CTCCCCATACAACTCAAAGTCCCCTGGGCCCAATTCAGAGTTATGTTTTTTTGGCACATACTAGAA
AGGCAGTGCCTCAGCCCTCCCTGAATCCATGGAGGTGTTCTGTTTGGGCTTTTACTGCTGCTGC
TCAGCTGGTTGCTTGAAGTACAGTAGGCCAGCCTGTTCTCTGCCATTCCCTAGTCATCCTGTGCCTCA
CCACAGCTTGCTTAGAGCAAGCCTTTTCTCAGACCTTAGGCACAGCCTCTCCTCTTTACCTGATCAATG
TTAAATGTAAGCACCCCTGATCCCAGGACATAAGGAAAGATGCCCAATTGTACTTTTGTCTATAGCCT
GTGAAATGGCTAGTTGATCATTTTTCCACAAAGAATTAGGTGTTAAGAGTTTTCTTCAGGCTTTACTT
AGGAGAATGGACTAAGCTGAAGGTGACTTCACCAGCAAGAGTCAACTCTAGAATTCAGGATGTTCTT
CTATTGTTTTCTTATCCATCTGTCAGGAAATGTAAGTTTGGTTTTATTTTGGCTTATCCAAGGGTA
AGCCAGAAAATAGAAATGATTATTTCTGATTAATAGCAGAAACTTTTCAATCTCAAATATATAAGGTG
TCTGCTCTTTTAAAGCTCTAAGCTAAGTCAAGAGCTAGGAAGCTGTTGATACAAATAAAAGTTTTTGAA
GGGA
ACGCGTAAGCGGCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018095.6</u>
Locus ID:	55709
MW:	28.6