

Product datasheet for **SC209451**

ARHGAP17 (NM_001006634) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ARHGAP17 (NM_001006634) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ARHGAP17
Synonyms: MST066; MST110; MSTP038; MSTP066; MSTP110; NADRIN; PP367; PP4534; RICH-1; RICH1; WBP15
ACCN: NM_001006634
Insert Size: 771 bp
Insert Sequence: >SC209451 3'UTR clone of NM_001006634
The sequence shown below is from the reference sequence of NM_001006634. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GACAATGATACCGAGAGCACTGCCCTGTGAAGAAAGCCCTTTCCAGCCCTCCACCACTTCCACCCTGG
CGAGTGGAGCAGGGGCAGGCGAACCTCTTCTTTCAGACCGAACAGTGAAAAGCTTTTCAGTGGAGGAC
AAAGGAGGGCCTCACTGTGCGGGACCTGGCCTTCTGCACGGCCCAAGGAGAACCTGGAGGCCACCACTA
AAGCTGAATGACCTGTGCTTGAAGAAGTTGGCTTCTTTACATGGGAAGGAAATCATGCCAAAAAAT
CCAAAACAAAGAAGTACCTGGAGTGGAGAGAGTATTCTGCTGAAACGCGCATAGGAAGCTTTTGTCCC
TGCTGTTAATGCGGGCAGCACCTACAGCAACTGGAATGAGTAAGAAGCAGTGCCTTAACATATCTATTT
AATAAAATGCGCTCATTATGCAAGTCGCCTACTCTCTGCTACCTGGACGTTTATTCTTATGTATTAGGA
GGGAGGCTGCGCTCCTTCAGACTTGCTGCAGAATCATTTTGTATCATGTATGGTCTGTGTCTCCCCAGT
CCCCTCAGAACCATGCCATGGATGGTACTGCTGGCTCTGTACCTCATCAAAGTGGATGTGACCCAT
GCCGCTCGTTGGATTGTCGGAATGTAGACAGAAATGACTGTTCTTTTTTTTTTTTTTAAACAATGTA
ATTGCTACTTGATAAGGACCGAACATTATTCTAGTTTCATGTTTAAATTTGAATTAATATATTCTGTGG
TTTATATGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001006634.3
Summary:	RICH1 is a GTPase-activating protein (GAP). GAPs stimulate the intrinsic GTP hydrolysis of small G proteins, such as RHOA (MIM 165390), RAC1 (MIM 602048), and CDC42 (MIM 116952). [supplied by OMIM, Apr 2004]
Locus ID:	55114
MW:	29.1