

Product datasheet for **SC209477**

Syntabulin (SYBU) (NM_001099747) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Syntabulin (SYBU) (NM_001099747) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SYBU
Synonyms:	GOLSYN; OCSYN; SNPHL
ACCN:	NM_001099747
Insert Size:	958 bp
Insert Sequence:	>SC209477 3'UTR clone of NM_001099747

The sequence shown below is from the reference sequence of NM_001099747. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CGCCGCACCGCCTTCCGTATCAAAACCTAATAGAAAGTTGTTGTTACCGTGTGCCAATGTGTCCCATGT
GGGTTGTGCCAGGTAGAGAAACAGGAAGTCAATCATCTGTGACAGTCTCTATTCTGTCGTTTTGCTCCT
TGGTATTTGATTTGCACTATATTTAGTTGAAGCCTGTTCACTGTTTAAAACCGGAGGTATCTTCAAAGG
CATGGAGACCTGGTTCAGTAAATGTCCACCAGTGGGGTATAGAAAGCATGCTCATGACCTGCGGTG
TCGTCTGAGGTACCCGTTCTTATCCTAGTGGTTCAGGAAGAGAAAACGCAGTTTGCACTTTCAAGACAG
CTTCTCTAAGGCTGGCATGTTATCTCCTTGCTTTGCTTTTTGCCGTTTTAAAATGTGTAATTGTTCCAG
CATTCCAATGGTCTTGTGCATAGCAGGGGACTGTAACCAAAAATAAACATGTATTTGTGTAATTGGTTT
GAAGAAGTCTTGAATAGCTCTTTACTGTCTTACTGGGGTTGATAAGATTTGAGTGTGCAATTTTTT
ACTAAATGTAGTCCAAAGTCTTAAATGGCTTGTGTTCTTAACTGTTAATTGATGAACTGTGCAT
AAGTTTACAATGTACTAATTATTTGCTTATTATATATAGTGTGTTTATTGGAAATTGTAACCACACAC
TTCAGCATGATGAAAATAAGATTAGTGTGTTCCATTTAAATAAATGTTTTATCCTCCCATAAAATAATA
ATTATTTCTGTGATTTTACTTCTATGTATAAGAAAGACAATGAATAGCAGATCTTACATCTTTATCAAT
AGGGGGTCTGTAAAAATTGAATTAATTAAGGAAGAGGGAAAAAGAGTTAGTTATTAAGAAAAATGAAGA
TTTCTACCCTTACTTAAAAAGCATATTTCTGAATAAAAGTAGATTTGTTTACTTTTTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001099747.2
Summary:	Syntabulin/GOLSYN is part of a kinesin motor-adaptor complex that is critical for the anterograde axonal transport of active zone components and contributes to activity-dependent presynaptic assembly during neuronal development (Cai et al., 2007 [PubMed 17611281]).[supplied by OMIM, Mar 2008]
Locus ID:	55638
MW:	36.9