

Product datasheet for **SC209505**

NDUFA4L2 (NM_020142) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: NDUFA4L2 (NM_020142) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: NDUFA4L2
Synonyms: COXFA4L2; NUOMS
ACCN: NM_020142
Insert Size: 690 bp
Insert Sequence: >SC209505 3'UTR clone of NM_020142
The sequence shown below is from the reference sequence of NM_020142. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCTGAAGAAGGACCGCCAGACTTCTAAGCCAGGCTGGGCTGCCAGTGCCATGCAAGCCACAGCCAG
CCAGCCCATCCACTTCTCCACTCCTCCCGCAGGCCCAAGGCATCACTCCGGCCACCCTGTCCCGT
ACTGCTTACACAGGCCGGTTCACAGCAGAGGGGAGGCTGCTCCACCCCTACTCTCCTCCCTTGCTCC
CAGCAGCGGAAGCGCCTCTGACCCTTGCTTGAAGTCCACGTGGGGGAGGAGGAGGCAGGCAGCACCAG
CAGGGGTCCACCAAGAGCCAGACCAGCCCTCTGCCCTCTACCCGGGCTCGAAGGGTGTGGCACAG
GCTACGTGTTGAGCGTGGCCTACGTGAGCCAACAAGAAGCAGGGGCTCTGAGTGCCAAGCGACGTGGC
GGGCTCCACGTTAGCCAGGCTCTGAGAGCCAGCCAGGGGCGGCGCTGCTCAGCTTGGGCTGGTCCAG
GGCCTGCCAGGCTGGGGCACCTTTGCCTCTGAGGCGCAGCGCACTCCTCCCCTGCCAAGCCTACTG
CCTCCCGCTGCCCGCAGTACCCCTCCAGCCCCACACCTGGGCCCTCCCCTGCCACTCCCCTCCCTTGC
TCCCCTCTGTCCCAGGGATCAAACAGAAGCAGCCGTGGGCAAAATACAATTTCAATTAACAAATTGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: NM_020142.4
Locus ID: 56901
MW: 24.1