

Product datasheet for SC209528

TOM1L1 (NM_005486) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TOM1L1 (NM_005486) Human 3' UTR Clone
Symbol:	TOM1L1
Synonyms:	SRCASM
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005486
Insert Size:	753 bp
Insert Sequence:	>SC209528 3'UTR clone of NM_005486 The sequence shown below is from the reference sequence of NM_005486. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CACAAAGGAGCTCAAAATGATGGTGACTGAGAAGAAAGTGGATGATCAGCTCACTACCACATCAAAGGT
GCCAACTCTCTAAACGTAGACTCTGTGCAGCTTTGAAGCCTGGAAGACAATACCTACCAACATGTCAA
AGCCATGGTGGCACATTTCTGCTATAATGAAGATTAATAGAATAACAGTTCCAGGATAACACTGATTC
CTGACAACAGCGTGAGATTTCAACAGAACTTGTTTGAACAAATACTCACTTAAACTTCAGCAGAAGA
AAAATTACTTAGTCCTTAGGCCAACCAATTTAACTGCAGTGTCATGTTTCACAGGCCTTCCTACATTTA
GAAATCGTCACACAGCTGTGATAAGAGTAGATTATTTACTATGAAATAATTCTGAATAGATGAAAGCA
TAAATGTGAGAACTGAATGTATTATTCAGGAAGAATACTGAGTGCCTTCATTTAACTAAAGTTGAAT
GTAAAAGTCAATTTGCACTTCTTTATAATCCTCTGGTTTAGAATTATAAATTGTTAAACCTTGATAAT
TGTCATTTAATTATTTTCAGGTGTCTGAACAGGTCCTAGACTCTACATTGGGCAGCCTTTAAATAT
GATTCCTTTGTAATGCTAAATAGCCTTTTTTCTCTTTTACTGCAACTTAATATTTCTATTTAGAACAC
AGAAAAAGAAAAATTTTGAATAAGTTGTACATTTGATGACAAATAAATCACTATTTAAACAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAAAGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_005486.3
Summary:	Probable adapter protein involved in signaling pathways. Interacts with the SH2 and SH3 domains of various signaling proteins when it is phosphorylated. May promote FYN activation, possibly by disrupting intramolecular SH3-dependent interactions (By similarity). [UniProtKB/Swiss-Prot Function]
Locus ID:	10040
MW:	29