

## Product datasheet for **SC209541**

### TPM1 (NM\_001018004) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** TPM1 (NM\_001018004) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** TPM1  
**Synonyms:** C15orf13; CMD1Y; CMH3; HEL-S-265; HTM-alpha; LVNC9; TMSA  
**ACCN:** NM\_001018004  
**Insert Size:** 773 bp  
**Insert Sequence:** >SC209541 3'UTR clone of NM\_001018004  
The sequence shown below is from the reference sequence of NM\_001018004. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGACTTTACTGGAGTTAAACAACATGTGAAAACCTCCTTAGCTGCGACCACATTCTTCGTTTTGTTT
TGTTTTGTTTTTAAACACCTGCTTACCCTTAAATGCAATTTATTTACTTTTACCCTGTCACAGA AAC
ATCCACAAGATACCAGCTAGGTCAGGGGTGGGAAAACACATACAAAAGGCAAGCCCATGTCAGGGC
GATCCTGGTTCAAATGTGCCATTTCCCGGGTTGATGCTGCCACACTTTGTAGAGAGTTTAGCAACACAG
TGTGCTTAGTCAGCGTAGGAATCCTCACTAAAGCAGAAGAAGTCCATTCAAAGTCCAATGATAGAGT
CAACAGGAAGGTTAATGTTGAAACACAATCAGGTGTGGATTGGTGTACTTTGAACAAAAGGTCCCCC
TGTGGTCTTTTGTCAACATTGTACAATGTAGAACTCTGTCCAACACTAATTTATTTGTCTTGAGTTT
TACTACAAGATGAGACTATGGATCCCGCATGCCTGAATCACTAAAGCCAAGGGTCTGTAAGCCACGCT
GCTCTTCCGAGACTTCCATTCTTTCTGATTGGCACACGTGCAGCTCATGACAATCTGTAGGATAACAA
TCAGTGTGGATTTCCACTCTTTTCAGTCTTCATGTTAAAGATTTAGACACCACATACAACCTGGTAAAG
GACGTTTTCTTGAGAGTTTTAACTATATGTAAACATTGTATAATGATATGGAATAAATGCACATTGTA
GGACATTTTCTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-Mlul  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



|                    |  |
|--------------------|--|
| <b>Components:</b> | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.   |
| <b>RefSeq:</b>     | <a href="#">NM_001018004.2</a>   |
| <b>Summary:</b>    | This gene is a member of the tropomyosin family of highly conserved, widely distributed actin-binding proteins involved in the contractile system of striated and smooth muscles and the cytoskeleton of non-muscle cells. Tropomyosin is composed of two alpha-helical chains arranged as a coiled-coil. It is polymerized end to end along the two grooves of actin filaments and provides stability to the filaments. The encoded protein is one type of alpha helical chain that forms the predominant tropomyosin of striated muscle, where it also functions in association with the troponin complex to regulate the calcium-dependent interaction of actin and myosin during muscle contraction. In smooth muscle and non-muscle cells, alternatively spliced transcript variants encoding a range of isoforms have been described. Mutations in this gene are associated with type 3 familial hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008] |
| <b>Locus ID:</b>   | 7168   |
| <b>MW:</b>         | 29.6   |