

Product datasheet for **SC209577**

FBXL14 (NM_152441) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBXL14 (NM_152441) Human 3' UTR Clone
Symbol:	FBXL14
Synonyms:	Fbl14
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152441
Insert Size:	785 bp
Insert Sequence:	>SC209577 3'UTR clone of NM_152441 The sequence shown below is from the reference sequence of NM_152441. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTGAGAACTCGGGGAAGCTCCAGAAGGTGAGGGAGAGGGGACAACGACATGGTTCCCGTGGATCTTTAA
CTTCCAGACTTGCCCGCTCTGCGCCTCCGGCACTCTGGTGATGACAGCTCAGGTTTCCCTGCCTGTCAC
TGCTCGGGCAGAGGCTGCTGCCAGGGCTTCTGCTCCGGTACCTTGGAAGCTGCATTCTCCTGCCGGT
TTCTCCAGTTCTGGGGACAGTGGTTTGCTCTGAGACCTCGCTTCTTTATGGATCCAAGGAGACTTGCT
TTTTCACTCTGTTCACTTTTTACTTGCTAGGATGGAATTGCAATTTGCAAGCTTCTTCGACAGGAAAC
TACAAGTTCCACACTTTAATTTTATACATATAAATATATACATGTGTACATATATCTATGTACAGGGGT
ATTATATATATACATATAAGATGATGATATATATAATGATGATATGTATTACTGAGAACGTAAAATATC
ATTACATAGTGATAGCTGGACACACAAGGAATTCACAACCCCCAAAGAAAATACATCTGGATGACCTG
CCTAGCAGTTTCCCCATGAGATAGAGGAATGTCTACGTATTTTCCTGTTTCTGCCCAGAAACAT
TTCAATCACTGACAAATCATTATCATTCAATAATGTTTACTGAGTGCCCATATGTGAAAGAAATCC
ACTCTACATTCACAGATGCATTTCTCTCCCCACGGGGTTTCCATTTAATGGGAACAATGTAGAATA
TATCTGTCTTCCCTTTATAGATTCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152441.3</u>
Summary:	Members of the F-box protein family, such as FBXL14, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]
Locus ID:	144699
MW:	29.8