

## Product datasheet for SC209600

### FAM40A (STRIP1) (NM\_033088) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | FAM40A (STRIP1) (NM_033088) Human 3' UTR Clone                                                                                                                                                                                                     |
| Symbol:                   | FAM40A                                                                                                                                                                                                                                             |
| Synonyms:                 | FAM40A; FAR11A                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_033088                                                                                                                                                                                                                                          |
| Insert Size:              | 759 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC209600 3'UTR clone of NM_033088</p> <p>The sequence shown below is from the reference sequence of NM_033088. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CCCATTTCCTGGGAAGAGCTGCTGCAGTAGGCTGTTGGTTAGGGGACTGAAATGGAGAGAAAAGATGA
TCTGAAGGTACCTGTGGGACTGCTCCTAGTTCATTGCTGCAGTGTCCCATCCCCACCAGGTGGCAGCA
CAGCCCCACTGTGTCTCCGAGTCTGCTCTGGGCTTGGGTGAGCCAGCTTGACCTCCCCTTGGTTCC
CAGGGTCTGCTCCGAAGCAGTCATCTGCTGAGATCCATTCTCCTTTACTTCCCCACCCTCCTC
TCTTGATATGGTTGGTTTGGCTCATTTCAATCAGCCCAAGGCTGGGAAAGCTGGAATGGGATGGG
AACCCCTCCGCCGTGCATCTGAATTCAGGGGTCTGCTGATGCCTCTCGAGACATACAAATCCTTGCT
TTGTGAGCTTGCAAAGGAGGAGAGTTTAGGATTAGGGCCAGGGCCAGAAAGTCGGTATCTTGTTGTGC
TCTGGGTGGGGTGGGGTGGTTCTGATGTTATCCAGCCTCCTGCTACATTATATCCAGAAGTAATTG
CGGAGGCTCCTTCAGTGCCTCAGCACTTTGATTTGGACAGGGACAAGGTAGGAAGAGAAGCTTCCCT
TAACCAGAGGGGCCATTTTCTTTGGCTTTCGAGGGCCTGTAAATATCTATATATAATTCTGTGTGT
ATTCTGTGTCATGTTGGGGTTTTAATGTGATTGTGATTCTGTTTACATTAAAAAGAAGCAAAATAA
ACGCGTAAGCGGCCCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

|                    |                                                                                                                                                                                                                                                                                           |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | SgfI-MluI                                                                                                                                                                                                                                                                                 |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |


[View online »](#)

|                    |                                                                                                                                                                                                                                                                                                         |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b> | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                                          |
| <b>RefSeq:</b>     | <u><a href="#">NM_033088.4</a></u>                                                                                                                                                                                                                                                                      |
| <b>Summary:</b>    | This gene encodes a member of the striatin-interacting phosphatase and kinase complex, which is involved in localization of the Golgi body. The encoded protein participates in cytoskeletal organization. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012] |
| <b>Locus ID:</b>   | 85369                                                                                                                                                                                                                                                                                                   |
| <b>MW:</b>         | 28                                                                                                                                                                                                                                                                                                      |