

## Product datasheet for **SC209602**

### **C1orf43 (NM\_001098616) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** C1orf43 (NM\_001098616) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** C1orf43  
**Synonyms:** HSPC012; NICE-3; NICE3; NS5ATP4; S863-3  
**ACCN:** NM\_001098616  
**Insert Size:** 763 bp  
**Insert Sequence:** >SC209602 3'UTR clone of NM\_001098616  
The sequence shown below is from the reference sequence of NM\_001098616. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACTATAACACATTGGAGAGTACTCTGTGACGGAGCTGAAGGACTCTTGCCGTAGATTAAGCCAGTCAG
TTGCAATGTGCAAGACAGGCTGCTTGCCGGCCGCCCTCGGAACATCTGGCCAGCAGGCCAGACTGT
ATCCATCCAAGTTCCCGTTGTATCCAGAGTTCTTAGAGCTTGTGTCTAAAGGGTAATTCCCAACCCCTT
CCTTATGAGCATTTTTAGAACATTGGCTAAGACTATTTTCCCCCAGTAGCGCTTTTTTCTGGATTGCA
TTCAGGTGTTATTCTTAATGTTTCTGTCAAAGCTTCTTAAAAATCTTCACTTGGTTTCAGCCATAGTTC
ACCTTCCCTGTTCCAGGTTTATTTAATCCAAAGGTGAGAGTTGGAGTGAGATGTCTCCATATCTATA
CCTTTGTGCACAGTTGAATGGGAACTGTTGGGTTTAGGGCATCTTAGAGTTGATTGATGGAAAAAGCA
GACAGGAACTGGTGGGAGGTCAAGTGGGAAAGTTGGTGAATGTGGAATAACTTACCTTTGTGCTCCACT
TAAACCAGATGTGTTGCAGCTTTCCTGACATGCAAGGATCTACTTAAATCCACACTCTCATTAAATAAA
TTGAATAAAAGGGAATGTTTTGGCACCTGATATAATCTGCCAGGCTATGTGACAGTAGGAAGGAATGGT
TTCCCCTAACAAAGCCAATGCACTGGTCTGACTTTATAAATTATTTAATAAAATGAACTATTATCAAAT
AAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001098616.3</a>
<b>Summary:</b>	General regulator of phagocytosis. Required to uptake Gram negative bacterium by macrophages.[UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	25912
<b>MW:</b>	28.6