

Product datasheet for SC209671

ZNF45 (NM_003425) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF45 (NM_003425) Human 3' UTR Clone
Symbol:	ZNF45
Synonyms:	KOX5; ZNF13
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003425
Insert Size:	793 bp
Insert Sequence:	>SC209671 3'UTR clone of NM_003425 The sequence shown below is from the reference sequence of NM_003425. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCAGAGGATTCACACAGGAAAACTCGATAAATATGTTTTACTATCTCAGATGGGTGCTGAAATATTTT
AATAATCAGAGCTATCATAGACAAAACATTTGTTTTATAGAGTCAGTAGTTCAGCCAAGTGATTGGGAG
ACCACACAGCAGAGAAGCCTCACAAGAGTGAGACATATGGACTGCATTGAGAACATTGACCATTAGCT
GATACATGCAGACAAGAGGATCAGGAAGGATGAGTCTGATCTGGAGTAAATCAGAAGTACTAAGATAGA
AATGCTGAATTCTGTTCCACTAGAATATAAGATCCAAGAGGGCAGGGACTTTGTTGACTGCCAAATCTA
CTCTGCCTTTTCAGTGCCTAATACGTTGTAATTTTTCAGTAGTATTTGAAATTACTGTGATATTTGAAA
TTCAGTAATATTTGACATTTTATTTATCTCTAGAAGTTTCTCTAAAATTGTACTCAGAAGAATTCTGC
AAGGCTTGGAGGATATATAAGTTAGTCATATGGCCTGATTTTCCATTTTGCAGATCTCGTGGACAAG
TGTTTATCAAAGTGAAGTTGCAGCTTGTTAGTGGGTTTCAGAAATCAGTTTCTGGCTAGCAACTAGAAG
TTTTGTATGTTAACAGTACTTTATTGAGGTGTAATCTGCATGCAGTAACATGCACAAGTCTCAAGTTCA
TGGCTTTCTGAATTTTGACAAATGTAATCATCACCAAGATCAATGTATACAATATGTCTATTAAACCAA
AAAGTTTTCTTGTGCCCATTTCAGTCAATCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_003425.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7596
MW:	31.1