

Product datasheet for **SC209747**

ZNF671 (NM_024833) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF671 (NM_024833) Human 3' UTR Clone
Symbol:	ZNF671
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024833
Insert Size:	760 bp
Insert Sequence:	>SC209747 3'UTR clone of NM_024833 The sequence shown below is from the reference sequence of NM_024833. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CAGAGGGTTCATGCTGGAGAAAAGCTTTAAGAGTGTAGCAAATGTGGGGGAAAGTCTTAGGCCAATGCC
CCTGACTTACTATATGGTGGGGAAGTAGCAGTAGTTAATGAGTGCAGCAGATGCAGGAAAGCCTTCCCC
TGGAGGCTGAACCTTACCCGCCATTGGGAATTTACACCGGACACAGGCCTTAGCAGTCTAAGCAATGT
GCTGTCTCTGTTCAAGCCCAACAGCTCACCTAGAGTGGAAGTCTGGGAGCAGCCATTGGGAGGGAACCA
TCAGTAAGAAGTGAACTTCATAGATATGGACATTCCCACTGGGGAGATTCCCTGTGAGTGTCAAGTAT
GTGAGATGCTTTCAGCAGCTGTGTTGCACTTTTAAATGGCTATTGGCCTTTGCTGGGGCAGGAGCCAT
CTGCTCTACCATCTGGCAGAATCATACTGCGTTTACCATTTACCCAGCATGCTTGTGACGGGCAGAC
CTCTCTTCTCTCCCCAGTCCCTAAAAGGTGTTGTGAGTGGTCTCACAGCCCACTAGGGGTCTTAATTTT
CTCTCTTTTGATGTAATGGCATGGAAATAATCAGCTTTGTTCAAGAGGACACAGAAGGATTCTGCAAA
TAGCCTGCAGAGACTTACCTGTGTTGATTGATTTTCATATGATGCTCGTTATGGATATATCCAATATCCA
AGTCAACCCAGCTCTGGAAGTGCCTGCTTACATTGCTCATGATAATAAAGGCCTGTTCTTTAAGCCACCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024833.3</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	79891
MW:	28.1