

Product datasheet for SC209805

KIAA0649 (PPP1R26) (NM_014811) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KIAA0649 (PPP1R26) (NM_014811) Human 3' UTR Clone
Symbol:	KIAA0649
Synonyms:	KIAA0649; NRBE3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014811
Insert Size:	782 bp
Insert Sequence:	>SC209805 3'UTR clone of NM_014811 The sequence shown below is from the reference sequence of NM_014811. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGTGGCGGCAGCTCAGTAGTGAAGGTCAGCCCTCGAGCTGTGGGTTTCGCGTCTGGGTTGCGTGCAT
TCGTGGAAAGCGGCGTAGCCGTGCGTGTGTGATGGTTCCGTGGCTGCAAGGAAGGGAACAGTCTAT
TATACATAGCCCTGTATATATGTACACCAACATGAACTTTTATTTAACAGACGTGTCCTGGTAAATAT
GATTTTTGTAGCTTTTTGTAAATTATTTAAAGTGATGTAAAAGATATTTTGGAAAATACTGTTGTTCA
ATTTTGTAGGGTGTTCCTAACTGCAGTTTTCTGTGTTCTGCATACAAGTCTTAGCTTAGGAAACATTG
GTTCTTACCATCACAGCCAGGTTACAGAGGCTCTGATGCTTTTGGTTGATTGCTGGTGAGAATGGCC
GGCGGTGGCTGCAGTGGTAGCCTGAGGAAGGCCAAGCTGCCCTCCCTGGGAATCACTCAGATGCCCAA
GATGTCCGTTGGGAAGCTCCCAGGACAGCACTTTTATACAGAGGACACCGCCTCGGCCCCACGTCCTT
AGAGGCCAGAGCACATCTGAAACTGCAATCACAGCCGTCCGCTGGAAAAACGTTTCGAAGCACAGTGG
CCAGCGAGCGAGCACGTGGCTACTCCCTGTTGCATGTCAAATCCACACAGATGTCAGCGGCAGTGCTC
GGCAGCCCAGCCCTGGGCTGGGTGGGTTTCATGTCCAATCTTGTTCGATCACTAGATGATTCTAACATC
GAAATAAACCTCTTTTATATGG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014811.5</u>
Summary:	Inhibits phosphatase activity of protein phosphatase 1 (PP1) complexes. May positively regulate cell proliferation.[UniProtKB/Swiss-Prot Function]
Locus ID:	9858
MW:	29.1