

## Product datasheet for **SC209839**

### CD168 (HMMR) (NM\_012485) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CD168 (HMMR) (NM_012485) Human 3' UTR Clone
Symbol:	CD168
Synonyms:	CD168; IHABP; RHAMM
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_012485
Insert Size:	807 bp
Insert Sequence:	>SC209839 3'UTR clone of NM_012485 The sequence shown below is from the reference sequence of NM_012485. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CCTATGGAGTGTCAAGAATCATGGAAGTAAACATCTGAGAAACCTGTTGAAGATTATTTTCATTTCGTCTT
GTTGTTATTGATGTTGCTGTTATTATATTTGACATGGGTATTTTATAATGTTGTATTTAATTTAACTG
CCAATCCTTAAATATGTGAAAGGAACATTTTTTACCAAAGTGTCTTTTGACATTTTATTTTTCTTGCA
AATACCTCCTCCCTAATGCTCACCTTTATCACCTCATTCTGAACCCTTTGCTGGCTTTCCAGCTTAGA
ATGCATCTCATCAACTTAAAAGTCAGTATCATATTATTATCCTCCTGTTCTGAAACCTTAGTTTCAAGA
GTCTAAACCCAGATTCTTCAGCTTGATCCTGGAGGTCTTTTCTAGTCTGAGCTTCTTTAGCTAGGCTA
AAACACCTTGGCTTGTTATTGCCTCTACTTTGATTCTGATAATGCTCACTTGGTCTACCTATTATCCT
TCTACTTGTCCAGTTCAAATAAGAAATAAGGACAAGCCTAACTTCATAGAAACCTCTCTATTTTTAATC
AGTTGTTAATAATTTACAGGTTCTTAGGCTCCATCCTGTTTGTATGAAATTATAATCTGTGGATTGGC
CTTTAAGCCTGCATTCTTAACAACTCTTCAGTTAATTTCTTAGATACACTAAAAAATCTGAGAACTCTA
CATGTAACATTTCTTCAGAGTTTGTATATACTGCTTGTCTGATCTGATGCTACTCAGCATTGATTA
ACATTTGTGTAATATGAAATAAAATTACACAGTAAGTCATTTAACCAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_012485.3</a></u>
<b>Summary:</b>	The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Dec 2008]
<b>Locus ID:</b>	3161
<b>MW:</b>	30.7