

Product datasheet for **SC209842**

CD168 (HMMR) (NM_012484) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD168 (HMMR) (NM_012484) Human 3' UTR Clone
Symbol:	CD168
Synonyms:	CD168; IHABP; RHAMM
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_012484
Insert Size:	807 bp
Insert Sequence:	>SC209842 3'UTR clone of NM_012484 The sequence shown below is from the reference sequence of NM_012484. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CCTATGGAGTGTCAAGAATCATGGAAGTAAACATCTGAGAAACCTGTTGAAGATTATTTTCATTTCGTCTT
GTTGTTATTGATGTTGCTGTTATTATATTTGACATGGGTATTTTATAATGTTGTATTTAATTTAACTG
CCAATCCTTAAATATGTGAAAGGAACATTTTTTACCAAAGTGTCTTTTGACATTTTATTTTTCTTGCA
AATACCTCCTCCCTAATGCTCACCTTTATCACCTCATTCTGAACCCTTTGCTGGCTTTCCAGCTTAGA
ATGCATCTCATCAACTTAAAAGTCAGTATCATATTATTATCCTCCTGTTCTGAAACCTTAGTTTCAAGA
GTCTAAACCCAGATTCTTCAGCTTGATCCTGGAGGTCTTTTCTAGTCTGAGCTTCTTTAGCTAGGCTA
AAACACCTTGGCTTGTTATTGCCTCTACTTTGATTCTGATAATGCTCACTTGGTCTACCTATTATCCT
TCTACTTGTCCAGTTCAAATAAGAAATAAGGACAAGCCTAACTTCATAGAAACCTCTCTATTTTTAATC
AGTTGTTAATAATTTACAGGTTCTTAGGCTCCATCCTGTTTGTATGAAATTATAATCTGTGGATTGGC
CTTTAAGCCTGCATTCTTAACAACTCTTCAGTTAATTTCTTAGATACACTAAAAAATCTGAGAACTCTA
CATGTAACATTTCTTCAGAGTTTGTATATACTGCTTGTATCTGCATGCTACTCAGCATTGATTA
ACATTTGTGTAATATGAAATAAAATTACACAGTAAGTCATTTAAACCAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_012484.3</u>
Summary:	The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Dec 2008]
Locus ID:	3161
MW:	30.7