

## Product datasheet for **SC209846**

### **NAT1 (NM\_001160174) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** NAT1 (NM\_001160174) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** NAT1  
**Synonyms:** AAC1; MNAT; NAT-1; NATI  
**ACCN:** NM\_001160174  
**Insert Size:** 799 bp  
**Insert Sequence:** >SC209846 3'UTR clone of NM\_001160174

The sequence shown below is from the reference sequence of NM\_001160174. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAACATGGTGATAGATTTTTACTATTAGATAAAGGAGTAAAACAATCTTGTCTATTTGTCATCCAGC
TCACCAGTTATCAACTGACGACCTATCATGTATCTTCTGTACCCTTACCTTATTTTGAAGAAAATCCTA
GACATCAAATCATTTACCTATAAAAAATGTCATCATATATAAATAAACAGCTTTTTAAAGAAACATAAC
CACAAACCTTTTCAAATAATAATAATAATAATAAAAAATGATTTTTAAAGATGGCCTGTGGTTATC
TTGGAAATTGGTGATTTATGCTAGAAAGCTTTTAAATGTTGGTTTATTGTTGAATTCCTAGAAAAGTTTT
ATTGGTAGATGAGTAAATAAAATATTGTAATAAACTTATTGTCTATAAAGTATATTAACATTGTTG
GCTAATAAATTTGAAAAAAGTGGTTTTTGGAAAGACTTAGGATATTATGGTGCTACATAATTTTTCC
TCGATGCTCTCTCCTCATCTTTCTGTCTCTAAATTAATTAATTAATTAATTAATTAATTAATTAATTA
AAGAATGAACATGAGCTTTTCTGTGTAGATCTGAGTTGAAATCCTGTGGACACTGGGCGAATTACTTT
TTAGATCTGTAGCTCTGACTCCTCAGGCATAAAATGGGAATAATGCTTTTACAGTTTAGTGCGGAACT
AAACTCCCAAATTTATTTGTTATATGGATCAAGTAATAACGTGAGTAATGTTTTGGTACAAAGTCATT
ATTTAATAAAAAGTTATTGTTCCATCTTGCTTGCCCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_001160174.3</a></u>
<b>Summary:</b>	This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
<b>Locus ID:</b>	9
<b>MW:</b>	31.3