

Product datasheet for **SC209847**

NAT1 (NM_001160173) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NAT1 (NM_001160173) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	NAT1
Synonyms:	AAC1; MNAT; NAT-1; NATI
ACCN:	NM_001160173
Insert Size:	799 bp
Insert Sequence:	>SC209847 3'UTR clone of NM_001160173

The sequence shown below is from the reference sequence of NM_001160173. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAACATGGTGATAGATTTTTACTATTAGATAAAGGAGTAAAACAATCTTGTCTATTTGTCATCCAGC
TCACCAGTTATCAACTGACGACCTATCATGTATCTTGTACCCTTACCTTATTTTGAAGAAAATCCTA
GACATCAAAATCATTTACCTATAAAAAATGTCATCATATATAAATAAACAGCTTTTTAAAGAAACATAAC
CACAAACCTTTTCAAATAATAATAATAATAATAAAAAATGATTTTTAAAGATGGCCTGTGGTTATC
TTGGAAATTGGTGATTTATGCTAGAAAGCTTTAATGTTGGTTTATTGTTGAATTCCTAGAAAAGTTTT
ATTGGTAGATGAGTAAATAAAATATTGTAATAAACTTATTGTCTATAAAGTATATTAACATTGTTG
GCTAATAAATTTGAAAAAAGTGGTTTTTGGAAAGACTTAGGATATTATGGTGCTACATAATTTTTCC
TCGATGCTCTCTCCTCATCTTTCTGTCTCTAAATTAATTAATTAATTAATTAATTAATTAATTAATTAAT
AAGAATGAACATGAGCTTTCTGTGTAGATCTGAGTTGAAATCCTGTGGACACTGGGCGAATTACTTT
TTAGATCTGTAGCTCTGACTCCTCAGGCATAAAATGGGAATAATGCTTTTACAGTTTAGTGCGGAACT
AAACTCCCAAATTTATTTGTTATATGGATCAAGTAATAACGTGAGTAATGTTTTGGTACAAAGTCATT
ATTTAATAAAAAGTTATTGTTCCATCTTGCTTGCCCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001160173.4
Summary:	This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Locus ID:	9
MW:	31.3