

Product datasheet for **SC209864**

CD21 (CR2) (NM_001006658) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD21 (CR2) (NM_001006658) Human 3' UTR Clone
Symbol:	CD21
Synonyms:	C3DR; CD21; CR; CVID7; SLEB9
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001006658
Insert Size:	799 bp
Insert Sequence:	>SC209864 3'UTR clone of NM_001006658 The sequence shown below is from the reference sequence of NM_001006658. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCTGTTGATCCATACAACCCAGCCAGCGATCAGAAGACAACTGGTGTGTGCCTCATTGCTTGAATT
CAGCGGAATATTGATTAGAAAGAACTGCTCTAATATCAGCAAGTCTTTATATGGCCTCAAGATCAA
TGAAATGATGTCATAAGCGATCACTTCCTATATGCACTTATTCTCAAGAAGAACATCTTTATGGTAAAG
ATGGGAGCCCAGTTTCACTGCCATATACTTTCAAGGACTTTCTGAAGCCTCACTTATGAGATGCCTGA
AGCCAGGCCATGGCTATAACAATTACATGGCTCTAAAAAGTTTTGCCCTTTTTAAGGAAGGCACTAAA
AAGAGCTGTCCTGGTATCTAGACCCATCTTCTTTTGAATCAGCATACTCAATGTTACTATCTGCTTT
TGGTTATAATGTGTTTTAATTCTAAAGTATGAAGCATTTTCTGGGGTTATGATGGCTTTACCTTTA
TTAGGAAGTATGGTTTTATTTGATAGTAGCTTCCTCCTCTGGTGGTGAATCATTTTCATTTTTACCC
TTACTTGGTTTGAGTTTCTCTCACATTACTGTATATACTTTGCCTTTCCATAATCACTCAGTGATTGCA
ATTTGCACAAGTTTTTTAAATTATGGGAATCAAGATTTAATCCTAGAGATTTGGTGTACAATTCAGGC
TTTGGATGTTTCTTTAGCAGTTTTGTGATAAGTTCTAGTTGCTTGTAATTTCACTTAATAATGTGTA
CATTAGTCATTCAATAAATTGTAATTGTAAGAAAAACATA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001006658.3</u>
Summary:	This gene encodes a membrane protein, which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Sep 2009]
Locus ID:	1380
MW:	30.2