

Product datasheet for **SC209877**

MCPIP1 (ZC3H12A) (NM_025079) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MCPIP1 (ZC3H12A) (NM_025079) Human 3' UTR Clone
Symbol:	MCPIP1
Synonyms:	dj423B22.1; MCPIP; MCPIP-1; MCPIP1; Reg1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_025079
Insert Size:	796 bp
Insert Sequence:	>SC209877 3'UTR clone of NM_025079 The sequence shown below is from the reference sequence of NM_025079. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCCTACAAGTCCCAGCACCCAGTGAGTAAGCTGCCTGTGGCTGGCAAGGGCAGCACCCAGCCTCCA
AGGGCCGTCAGGCTGGGCTTTGGGCCATTGAGCAGCCATTCCCAGCCCTGAGGCCACCCAGAGGCT
GGACAGAGGGAGGATTCAAGTCGGGAAGGAAACCCACAAACCAAGATACTGTAGGATTGGTTCTGGCC
CATGCAGCACCTCTAGCTGTCTGCCTCAGTGGGTGAGAAGCGATCACCTGTTGATACACATTGTATCT
CTGTAGTTTAAGGAGACGCTGCCGTAACGGCGTCGGTCCGTGGCTGAGGCCAAACCGTCTTTCTCT
CAGAGGGTGGGGAGGGAGGTGGGGCAGCAGAGGCCTGGGCTGGGTGCCCTGTGCACGCCACCCACTT
CCGCCCTACCCCTGGGACGTTGGCCTTGGCTGGCTAGTTGGGCACCGTGTGCCTGCCCTCCAAGGCCT
CCTCTACGCCAATGAGGCCTCATCTGTCTCTCGCTGGGCACGTGGCTTCATGTCAGTAAGCAAGATGC
TTCTTAATAACCCACCTTCTGCCCCACTCTATTCTTATCCTGCTGCCCTGTAGGGGTCAAGGGCCCT
CCGTCTACACCCTCTTCTCTCTCCATCCTTTATTCAGAGTCATCTCGCCCTTCCCCATGGGTGGGG
AACCTGTGTTTGTGTGTGCACATGTAATTTTAAATATTTTAAAGCAGAAAGTCCTTACCTCCTGTAA
CACATCAATAAAGTACAATCATTGTGAGCCCTTTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_025079.3</u>
Summary:	ZC3H12A is an MCP1 (CCL2; MIM 158105)-induced protein that acts as a transcriptional activator and causes cell death of cardiomyocytes, possibly via induction of genes associated with apoptosis.[supplied by OMIM, Mar 2008]
Locus ID:	80149
MW:	28.2