

Product datasheet for **SC209918**

PPM1M (NM_001122870) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PPM1M (NM_001122870) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PPM1M
Synonyms: PP2C-eta; PP2CE; PP2Ceta
ACCN: NM_001122870
Insert Size: 815 bp
Insert Sequence: >SC209918 3'UTR clone of NM_001122870
The sequence shown below is from the reference sequence of NM_001122870. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGTCAGGGCCAAGAGAGCAGTGACCACTGAGGATTGACACTGTATCCCAGAAGTCTCTAGTGCCCG
GGTGTGGTCTGGGCATCCCTCCAGTGTGACCAAGAGCAAATCCTGCCTGCCCTATCCCTAGCCACCGCC
CAGTGCTCTCACTATCCACCTCAACACACATCCATCTCAAGAGGAACATTTATACCAGGCAGTCAGAGC
TGGAAGTGTATGGAGAGCCAGCCACCAAGTCTGCCTTTTGGCGGTGATAACCTTCTCTGGCAGAGTG
ACTTTACAACCTAACTAGGAAACCCATGTGAGGCTCCTCAGACAGGATCTTGAACAGCCCAAAGTATCA
TTCTCAGATAGGGGCACCCAAGCTAAGGGTATTAGCCAAAGATGCCAGGATGGGTAGCTAGCCCATGTT
TAGATCCAGGTCTCCAATTCATGGTTATCAGGGCATGTGTTCAACAACCCCAAAGTCCACGCAGGTGG
CTTGTAGAAACCTTTGGGCAGCCTCATGTCTGCTAAAACAGCCATCTTCAAGACAGCCCTGAAAAGAG
ACCAGTTCAGGTCTGCCCTGCTGTTCTTTGCTGGAGATGAGGAACAGGTGCTGGGGCTAAAGTTTGGG
GTAGAGCACAAGGACAAGAGGAACTCTTGGAGTTGGCTGGGTGAGAGGGCTCTCCATTTGCTACCTGT
AGTAGCCTGCCTTTAACTGTTGCTTCTCCCTAGTTCAGCCCTGCCCTGGTCTGATGCCCAACACT
GCCCTTGCTTTGTTTTCCCTGTACCTCCCTATTATTAATGTTTTCTACAGAAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001122870.3](#)

Locus ID: 132160

MW: 29.8