

Product datasheet for SC209942

C8orf34 (NM_052958) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C8orf34 (NM_052958) Human 3' UTR Clone
Symbol:	C8orf34
Synonyms:	VEST-1; VEST1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052958
Insert Size:	807 bp
Insert Sequence:	>SC209942 3'UTR clone of NM_052958 The sequence shown below is from the reference sequence of NM_052958. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGTCTACGCTGGTCTACTCTGGTTTGTAACAGAGAGAAGAAGTTGCAAGTGGTCCTTAAAGAAATG
CAGTTATTCAAATCCTTATGTATAGTTCTAATTTTATTTACTATGTGAATCATTTAGAGAATGGTTAT
TTTTATAATCTCAATAATAAACAAGTACTATCGTAGCCAGGGGGTGCCCAAGTATGTAATCAGAGAACA
CTAATCCTGGCAAAGGATTGTGGGGTGGTCAGGAGGCCGCTGCCTTTTGACATGGTTAGAGTCCTGGG
TTTAGATTACTTTATAAATAGTTATTAAGATTAATATTTGATATATTAATATTGCCTGATTCAAATAA
CTTTGCTTAAATTTTTCTGTATTTAACTTGTCTAGGGTTTTCTACTTCTTCTATTTTTATAATATTT
TCTCATTTTTTAGATTACTAAACATTCTACTTATCCACATGTGCCAAGATGTGCCATTTAGGTCTCTGC
AGAAATTAGCTGTTTGGTATCATCACCTGCCATACATGTTGTGTGCAGCAGTGAGAAGTATATTCTCTA
CCACTGCTCTTTCTCATGATTATTTCTCTATTCTTATTAATGAATAATTGATATGTTAATTAAGTA
AAGCAATCAATAAGGAATGAAATGCATGGCCCAAAACATAGAAAATGAAAATCGTTTGACCATTTATC
CAATTATTGAACCAATTAAGATGATCCATTTTAAATTTAATATATGCTCGCTTTTGAATGTGTTATG
CATTGATGTTAAGTGAATGGATCTAAGTAAATCTATTGATATCTCTAC
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_052958.4
Summary:	This gene encodes a protein that is related to the cyclic AMP dependent protein kinase regulators. Naturally occurring mutations in this gene are associated with an increased risk for severe toxicities, such as diarrhea and neutropenia, in patients undergoing chemotherapeutic treatment. [provided by RefSeq, Mar 2017]
Locus ID:	116328
MW:	31.1