

Product datasheet for **SC209952**

PWWP domain containing protein 2B (PWWP2B) (NM_001098637) Human 3' UTR Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | 3' UTR Clones |
| Product Name: | PWWP domain containing protein 2B (PWWP2B) (NM_001098637) Human 3' UTR Clone |
| Symbol: | PWWP domain containing protein 2B |
| Synonyms: | bA432J24.1; pp8607; PWWP2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pMirTarget (PS100062) |
| ACCN: | NM_001098637 |
| Insert Size: | 815 bp |
| Insert Sequence: | >SC209952 3'UTR clone of NM_001098637 The sequence shown below is from the reference sequence of NM_001098637. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGGGACATCGTCTGGGGTCACCGACGATGAGGGCGCACCTGCTGCTCTGGAGCTTCCGATCTCTGTTTCG
GGGACCTGTGAGCCTTCTGGCCGGCTGCGTGCAGAGCCCACTGGGCACGGTGGTGGCCTGGTGTGAG
GCCCGGGGGACCGGCAGTGTGTCCAGGGAGGGGATGGCCCTGAGCCCAGCGGCTCCTCCCGCTGAGT
GTATTTCTTCCACCACTCAGCGCATGGCTGCCCTGGCTCACGAGGCAGTCGGGACTCGTGACTTGCTG
GCTGCTGGCTGCTGCTGCCTCGCGCGCTGGGGTCCATGGAGGAACTGGGCTGCCGCCCGGCCTCAC
CTGCTGCCCGCATGCTGGGGTCCCATGGAGGAACAGGCCTGGCGCCCGGCCTCGCCACGGGCTGGT
GTGGGCAGCTGTTCCCTGCCTCGCAGTGTCTGGAGGCAGCTGTCTCGCAGACTCAGCTTGGTCTCCCG
CAGGCTTCAGAAAAACCAATTGCACGTGTGGGATTCTTCCCAAGCTTCCAGAGCGTGGGGCCCGCTGG
ACTTGAGAGCTGGGGACCCACCCACCGTCTCCCACTCCAAGTTCATCCAGAGCGTGGGGCCCGCTGG
GCACCTGGCACCTGGCCTGCAGATGCTGCTACGAGTGACCAGTGTGTGGAGGAGGCCAAGGGCGGC
CCCTGGGAAACCGGGCTTCAACAGTACAAGGAAAAAGAACTTGTCTGTTTTGTAAGAGCGTGTTCCTT
TTCCTCTTTTTTGTCACTTTAAAGACCAAAAAGAATAAAGAAAAGAAAAAGAAAAA
ACGCGTAAGCGCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

| | |
|------------------------|--|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_001098637.2</u> |
| Locus ID: | 170394 |
| MW: | 29.2 |