

## Product datasheet for SC209954

### OSCP (ITGBL1) (NM\_004791) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | OSCP (ITGBL1) (NM_004791) Human 3' UTR Clone                                                                                                                                                                                                       |
| Symbol:                   | OSCP                                                                                                                                                                                                                                               |
| Synonyms:                 | OSCP; TIED                                                                                                                                                                                                                                         |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_004791                                                                                                                                                                                                                                          |
| Insert Size:              | 822 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC209954 3'UTR clone of NM_004791</p> <p>The sequence shown below is from the reference sequence of NM_004791. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAAATCTGGCTTGGCTCAGAATATCCTTAAACAATTACATGAGAGAGGTCTGGATTCTTATTTTTTCTGG
GCCATTAGAACAGATAAATGCGAAGGAAACCATGTATATTCACCACTAGGACAGGTAAAAAGACCATT
GTATGTTTTTCTATTTCTGAATTACGAATGAAATCCGAGTACCTATTAGAAATGAGTTATGCAAATTTA
GATGCAAATAACATTAGAAAAAAGATTCTTCCATAATTAACATAAGTGGTTCCTAACGAGAGCAATT
TTTCCACCCAAAAGTCATTTGGCAACATCTACAGACAATTTTGATTGTCACACTGGGTGGGTAGGAAG
GTATGCTGCAGACATTTGGTGGGTAGAGGCCAGGGATGCTGCTGAGCATCCCGCAGTGTACAGGACAGC
CCCCAAACAAGGAATTATCCAGCCCCAAATGCCAATAGGGCTCAAAGTGAAGAACATTGAGTTATATGG
CTATTAGAAATCCACATTCTTACACAAGAAAGACCATATTAGAATCTAAGGAAAACATGCATATTCACA
TTAATTAATCGATCAGATTTTCCAGAATTCGTATCAGTCACCATTTTAAATATGGGGACAATGAAGAC
AAGCACACAGGAGGTAGAATATCAGAGTGGGGCTGGATCAAGGGCAAAACTGGTCATTAAGTCATCTG
ACATTAATCATTAGCCACTAAGTTATTTGTCTACTCTCACTTTAAACTCACCAAGAAGATTCTCTT
AAAGAAATTATGAAAAATGTACAATTTAAATTTAAATAGTGACAGAAAGTTGTTTATA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

|                 |                                                                                                                                                                                                                                                                                            |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components:     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| RefSeq:         | <u><a href="#">NM_004791.3</a></u>                                                                                                                                                                                                                                                         |
| Summary:        | This gene encodes a beta integrin-related protein that is a member of the EGF-like protein family. The encoded protein contains integrin-like cysteine-rich repeats. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]                          |
| Locus ID:       | 9358                                                                                                                                                                                                                                                                                       |
| MW:             | 32.1                                                                                                                                                                                                                                                                                       |