

Product datasheet for **SC210013**

DATF1 (DIDO1) (NM_022105) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	DATF1
Synonyms:	BYE1; C20orf158; DATF-1; DATF1; DIDO2; DIDO3; DIO-1; DIO1; dJ885L7.8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_022105
Insert Size:	810 bp
Insert Sequence:	>SC210013 3'UTR clone of NM_022105 The sequence shown below is from the reference sequence of NM_022105. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCAGGACTGGGAGTTGCAGCTTTGTGTAAAGCTGATCACAGACACCGGCTGCACCATCAGCGGGAAGCA GAGCCCATGTCCAGGATGCCTCCTGCTGCCCTGTGTCCATCCCTAGTCTGTCAGGACTTCCTGTCACTG TTTTCCAAAGCTGTAAACCTCACTGGTGAACGTTACCTTAATGATTGATTCTTTAATCTCTGTTTTCA CTCTCAGGCTCTGGTAAGTATTTCGTATTCTCTTCATCCAGTCTGATTGCATAGCCACACTGCCCGGCA CGCCACATCCACCCTGTCTGCACATGAGTTGTTCTGACAACAGCGCTGTATACGCTTCAGTTTTTCCA CATTGTCCACGGCCAGCACATGAAAGCATCACTTCTTTTTATGTTGTGGGAATCTTTGCAAGTTAGTG TTGCATCTGATTTTCAGGTGTACATTTATTTTTGACTGGGCAGATAGGGGATTTTTTTTTTCCATGTC CGATTACACGCTACACACCCACATGAACACATTGAACTTCGAAGGCCACACACTCCTGCTTCATAGG CCCCACGGTAAGTGAGTTCACACCTAGAACACTGTCCTGACCGCAGGACGCGTGCCTTGGACTTGGTAT TCTACATGTGACTGGCTTTCTTGCCCTCGTCTCTTGAATGTTAGACTCTTAAGATCATATCCTGCCCC AAATTTCAAATTAATGAAATGAAGATATTTCAAACAGATCTTTGAAACCTCAGATTCTGTGGTGCAATT TTAATGTTTTCTTGTCTCAGTTTTCTGCTATAAACTATTTTCAATTCA AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCG CAACCTGCCATCACGAGATTTGATTCCACCGCCGC </pre>
Restriction Sites:	SgfI-RsrII



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_022105.5
Summary:	Apoptosis, a major form of cell death, is an efficient mechanism for eliminating unwanted cells and is of central importance for development and homeostasis in metazoan animals. In mice, the death inducer-oblierator-1 gene is upregulated by apoptotic signals and encodes a cytoplasmic protein that translocates to the nucleus upon apoptotic signal activation. When overexpressed, the mouse protein induced apoptosis in cell lines growing in vitro. This gene is similar to the mouse gene and therefore is thought to be involved in apoptosis. Alternatively spliced transcripts have been found for this gene, encoding multiple isoforms. [provided by RefSeq, Jul 2008]
Locus ID:	11083
MW:	30.7