

Product datasheet for **SC210016**

RSRC2 (NM_198261) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RSRC2
Synonyms:	arginine/serine-rich coiled-coil 2; FLJ11021
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_198261
Insert Size:	831 bp
Insert Sequence:	>SC210016 3'UTR clone of NM_198261 The sequence shown below is from the reference sequence of NM_198261. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TCTTCAATGCGAGGAATGGATGCAGTTGAAAAATGATCACACTTGTAAGTTTGGGACTTATAGACTTC TTGTTCTGATGTCACGTCCTTGTTCACCAACAGCTAGCACTCTAGCTTGCATGGGTGTTGCATTGACT TTAATTTATTGAAAAATACAAATTTTGTAAATATCAGATCAGTGATACTGGTGTAGTGTGTAATCA GGTTAAACCCACTTCCATTAACTTGACAGGACTATAGAAGGATAATATTTTTTAGTTCATGAATTCTA CTTTTCAAATATATAAAAGCTGCAGGTGGGGATAAAATCTCATACATGGATTTTTTCGTGTCGCTGTC TTGTGTACTTTTGTACTTAACCTTGTACAGTTATTTTCATCTCTTGAAACATGAAAGAAATGTTATGTA GATGTTCTTTAGAAGATCTGGCCATTTGGTACATAATCCAGCACAGATAAGCTGGGTGGTAATGATAAT AAAAATGGTTTTCTCAAACTGGTGTAAATTTAAGTTACCTGGGATGTTTCTTTGAATTTGTTTTATAG TTTCTGTAGCATTGGCAATTGCTGTTAGAAAACACTAGCTAGAAATCCCCTCCCCACCACCTTTTTTA AGGCCAGTTAACTATACTACAGTCAATACCGTGGTGAGCAAAATGTAAAGGTGGAAGGAGAAAACCT ACTAAATAGTATGTTTTCTATTATAAGGGACAGACTTGGTATTCAGTATTTGTCAAATATTACATGT GTTATTTCAGGAGATAGATTAATGCATTAAAGGGATGTAAGCACTTTTATTTTAATAAAGTGCCTTATAA CAA ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198261.2</u>
Locus ID:	65117
MW:	31.2