

Product datasheet for **SC210024**

WDR20 (NM_181291) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: WDR20

Synonyms: Bun107; DMR

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_181291

Insert Size: 820 bp

Insert Sequence: >SC210024 3'UTR clone of NM_181291
The sequence shown below is from the reference sequence of NM_181291. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCAGTTCTCCAGGTGGAAGTGTAGTGTAGCGACCTCACTGCTGCGCGCACAGTCTCCCGGACTTGA
CTCGAGGGAGTGACGAGGAGGAGCTCCGAGCTGCGCCTGAGCCGTGCCAGCCGGCGGACCTCAGGCGGT
GGACGTCGGCGATAGCCGTGTGGACGGTGACCGGCTCACTCTGCGGCCCGTGTCCCGTGTCTACCC
AAAGAAGTTGTTTCCATTTTAAACCGGTCTTTTGGGGCTGCAGTAAAAATAAGAAATGGAGTTTCTT
GCTTTTACTCTAAAAATCAATGTAATTAATTTATATATATAATATATACATATATACATAGTGT
AAAATAAAATGTTTCTTGACAAGAAATCCCTGAAATTCAGCTGTTATAGTGCTTCACTGTTTTGCA
CTGATTTTCTATACCTTAGGTGGTCAGAAGACAACCTGAATGCACTCATAGAGAAACTGTTACTTT
CTGACGTAATGTAATTCAGGAAGACAGACGCTGCAATCACAGATTTTAAAAAATGTTTGCACTTAAAA
ATAGTTGAATGCTGGTGAAAGTTACTTTGCAGATGGGTGTAAGGACTCATGGCCCTCTGAGGTGCGGC
GTGAAGATGCCCTTTTACCCGTTGACGTTATTTACGTAATAAATGTTGTTTCCAATGCAATCA
ACTCTGTATTATATGTATAAATATTGTAATTCGCAATTGGGAAAAATAGTTACTTCACTAGTAATTT
CATCATTTAAGAGTGATTTTCTAATTCACAAAAGTTAATATTAACACTATCTTGAATATA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_181291.3
Summary:	This gene encodes a WD repeat-containing protein that functions to preserve and regulate the activity of the USP12-UAF1 deubiquitinating enzyme complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jun 2011]
Locus ID:	91833
MW:	31.5