

Product datasheet for **SC210091**

THOC2 (NM_001081550) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: THOC2 (NM_001081550) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: THOC2
Synonyms: CXorf3; dj506G2.1; hTREX120; MRX12; MRX35; THO2
ACCN: NM_001081550
Insert Size: 818 bp
Insert Sequence: >SC210091 3'UTR clone of NM_001081550
The sequence shown below is from the reference sequence of NM_001081550. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CATCATAAGTCTCGGACAAGCACAGATAATGGAAGACTTCCATCAAGGTGAGATCGGACTGGAAGTGT
TCGGCTGCGACCAGAAATTTATTTCTGAGTAAATTGCCGAGAATTAAGAATGAAGAGGGCCATTGTC
ATCTCCTTAAATTATTCAGTTACCTGCTTTATTGCTCCATGTGGAAAACCTAAAATTGTTAAGTTGTGC
ATTACTGTATTTAACTTGTGCTTAGTTTCTACATGTTTATTTTCAGTAATGGCTGAAAGTGTAACT
GTTCCATACTTTTAGCACAATGTGCTGCATAAGGTTACCTGTGTACAGAGTTTACTTTAGATTAACATA
AATATTGCCTGGGTTTCAGTTTTTATTTCCATTCTGAAATGCTTCTTTTTATTGTTTGAAACTGAAAAT
AAACAATTGTTGAACCCTTTTGATTTTACCTCATTTTAAAAGTGTTTAATTTATTATTTGGCTTGTTC
TTAATATTAGTCACTAAAAGCAGTGGGAGCATTGTCTTATGAAATGCTTAGGAATCATTTTATATAGTA
CATGTACAACATTAACGTGTTTAAAAAAGAAAAAGGTACCAGCGATCACTTGCCCTTGCCATTTTTT
CTTGTAATTATGTTAGACAAATCTTGGCGCGGGGGGATCAAAACATAATTGTTTTAATTCTACAGCTG
TAGGAGCTTTGTATTGCTGAACTTTCATCTGAAAAGTTTCACAGTGACATTTTTAAAAGAGAATTTTT
TTATCTGCCGAATTCTACCAGTGTAACTTTTTTCTAAATAAAACAATAGTTTTTCTCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001081550.2
Summary:	The TREX multiprotein complex binds specifically to spliced mRNAs to facilitate mRNA export. The protein encoded by this gene is a member of the THO complex, a subset of the TREX complex. The encoded protein interacts with the THOC1 protein.[provided by RefSeq, Jun 2010]
Locus ID:	57187
MW:	31.6