

Product datasheet for **SC210102**

MS4A4A (NM_148975) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	MS4A4A
Synonyms:	4SPAN1; CD20-L1; CD20L1; HDCME31P; MS4A4; MS4A7
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_148975
Insert Size:	822 bp
Insert Sequence:	>SC210102 3'UTR clone of NM_148975 The sequence shown below is from the reference sequence of NM_148975. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAAATCAAGCGATCGCC GCATCTCCACACCACTTAATGAGGTTTGGGCCACCAAAAGATCAACAGACAAATGCTCCAGAAATCT ATGCTGACTGTGACACAAGAGCCTCACATGAGAAATTACAGTATCCAACCTCGATACTGATAGACTTG TTGATATTATTATTATATGTAATCCAATTATGAAGTGTGTGTATAGAGAGATAATAAATTCAAAATT ATGTTCTCATTTTTTCCCTGGAACTCAATAACTCATTTCACTGGCTCTTTATCGAGAGTACTAGAAGT TAAATTAATAAATAATGCATTTAATGAGGCAACAGCACTTGAAGTTTTTCATTCATCATAAGAACTTT ATATAAAGGCATTACATTGGCAAATAAGGTTTGAAGCAGAAGAGCAAAAAAAGATATTGTTAAATG AGGCCTCCATGCAAAACACATACTTCCCTCCCATTTATTTAACTTTTTTTCTCTACCTATGGGGAC CAAAGTGCTTTTTCTTCAGGAAGTGGAGATGCATGGCCATCTCCCTCCCTTTTTCTTCTCCTGCT TTTCTTTCCCATAGAAAGTACCTTGAAGTAGCACAGCCCGTCCTTGCATGTGCACGAGCTATCATTTG AGTAAAGTATACATGGAGTAAAAATCATATTAAGCATCAGATTCAACTTATATTTCTATTTTCATCTT CTTCTTTCCCTTCTCCACCTTCTACTGGGCATAATTATATCTTAATCATATATGGAATGTGCAACA TATGGTATTGTAAATACGTTTGTGTTTATTGCAGAGCAAAAAATAATCAAATTAGAAGCAA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_148975.3
Summary:	This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features, similar intron/exon splice boundaries, and display unique expression patterns in hematopoietic cells and nonlymphoid tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Locus ID:	51338
MW:	32.1