

## Product datasheet for **SC210130**

### STAG3L1 (NM\_018991) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | STAG3L1 (NM_018991) Human 3' UTR Clone                                                                                                                                                                                                             |
| Symbol:                   | STAG3L1                                                                                                                                                                                                                                            |
| Synonyms:                 | DKFZp434A0131; FLJ78546; MGC40269; MGC131759; STAG3L2; STAG3L3                                                                                                                                                                                     |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_018991                                                                                                                                                                                                                                          |
| Insert Size:              | 847 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC210130 3'UTR clone of NM_018991</p> <p>The sequence shown below is from the reference sequence of NM_018991. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GGGCCAGGGACCTTCCACCTCCTAGGGTAAACCAGGAGAGATTGCTTGCTTCACTTGTACAAGGCAGG
AACGGTGGCATGGGGTGGGGGAACTTGGAGTTGGAAGGTGGCTAATCTTTGATTCTATGTTTTGATC
CTCCTGGCACTCCAGACCTGGGTGATGTGCAGGAGAGCACACTGATAGAAATCCTTGTGTCCAGTGCCC
AGCAACTCCTGCCTCAGCCTCCCGAGCAGCTGGGACTACAGGCGCCGCCACCACGCCTGGCTAACTTT
TTTGATTTTTAGTAGAGACGGGTTTTACCGTGTTGGCCAGGATGGTCTTGATCTCTTGACCTTGTA
TCCACCTGCCTCATCATCCCAAAGTGCTGGGATTACAGGCGTGAGCCACTGCGCCAGCATGTTAGACA
ATTTTTAATTCATCCTCTCTGTGCTGTTGTTTTCTCAGCTGTGAAAGGAATATTCTGGTGGGGACAAGG
TTACAGAGTTGCTGAGAGGGTCTCATGACATGAAGGTAAGTGGCCTTGGCACAGTGCTGGGGGGCGGG
GACTCCGCACATGCCTGTGATGTACAGTTACTGTCAAGTTCACAGTGAACCTTCCTCTTTCTGTT
GACTTTCCCACTCCTGTAACCTCCCTCCCTCCCTCTCTCTCTCTCTCACTCACGCACAC
GCACACACACACACACACACACACACTCCATTCACTGTCTCCATGACTCTGGAGTAAACTAAC
GTCTCGAGTTGCCATTGGAAGCCCCATTGTCCTATTTGTTTTGTTTGTGTTGAGACGGGTCT
CGCTCTGTCCCCAGGCTG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

|                 |                                                                                                                                                                                                                                                                                            |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components:     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| RefSeq:         | <u>NM_018991.2</u>                                                                                                                                                                                                                                                                         |
| Locus ID:       | 54441                                                                                                                                                                                                                                                                                      |
| MW:             | 31.6                                                                                                                                                                                                                                                                                       |