

## Product datasheet for **SC210158**

### SERPINE2 (NM\_006216) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	SERPINE2 (NM_006216) Human 3' UTR Clone
Symbol:	SERPINE2
Synonyms:	GDN; GDNPF; PI-7; PI7; PN-1; PN1; PNI
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006216
Insert Size:	852 bp
Insert Sequence:	>SC210158 3'UTR clone of NM_006216 The sequence shown below is from the reference sequence of NM_006216. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTATTCATGGGGCAGATAAACAAACCCTGAAGAGTATACAAAAGAAACCATGCAAAGCAACGACTACTT
TGCTACGAAGAAAGACTCCTTTCCTGCATCTTTCATAGTTCTGTTAAATATTTTGTACATCGTCTCTT
TTTCAAACACTAGTTCTTAGGAACAGACTCGATGCAAGTGTCTGTTCTGGGAGGTATTGGAGGGAAAA
AACAAAGCAGGATGGCTGGAACACTGTACTGAGGAATGAATAGAAAGGCTCCAGATGTCTAAAAGATTC
TTTAAACTACTGAACTGTACCTAGGTTAAACAACCTGTTGAGTATTTGCTGTTTGTCCAGTTCAGGAA
TTTTTGTTTTGTCTATATGTGCGGCTTTTCAGAAGAAATTTAATCAGTGTGACAGAAAAAAA
TGTTTTATGGTAGCTTTTACTTTTTATGAAAAAAATTTTGCCTTTTAAATCTTTTCCCCATCC
CCCTCAAAGTCTTGATAGCAAGCGTTATTTGGTGGTAGAAACGGTGAAATCTCTAGCCTCTTTGTGT
TTTTGTTGTTGTTGTTGTTGTTTATATAATGCATGTATTCACTAAAATAAAATTTAAAAACTCC
TGTCTTGCTAGACAAGTTGCTGTTGTGCAGTGTGCCTGTCACTACTGGTCTGTACTCCTTGGATTTGC
ATTTTTGTATTTTGTACAAAGTAAAAATAAACTGTTATGAGTAGTAAAAATAAAGCTATTTCTCTGCTA
TTTGAAAATACAATAGAAGAACTGAGCCTTTTAGACATTCGTCAGCCTCTTCTAATAAACCTTTGTAC
TATGTAACATCAGGAAATTCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_006216.4</a></u>
<b>Summary:</b>	This gene encodes a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. Thrombin, urokinase, plasmin and trypsin are among the proteases that this family member can inhibit. This gene is a susceptibility gene for chronic obstructive pulmonary disease and for emphysema. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]
<b>Locus ID:</b>	5270
<b>MW:</b>	32.8