

## Product datasheet for **SC210160**

### REXO1 (NM\_020695) Human 3' UTR Clone

#### Product data:

|                |                                   |
|----------------|-----------------------------------|
| Product Type:  | 3' UTR Clones                     |
| Symbol:        | REXO1                             |
| Synonyms:      | EloA-BP1; ELOABP1; REX1; TCEB3BP1 |
| Mammalian Cell | Neomycin                          |
| Selection:     |                                   |
| Vector:        | pMirTarget (PSI00062)             |
| ACCN:          | NM_020695                         |
| Insert Size:   | 848 bp                            |

**Insert Sequence:** >SC210160 3'UTR clone of NM\_020695  
The sequence shown below is from the reference sequence of NM\_020695. The complete sequence of this clone may contain minor differences, such as SNPs.  
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTTTCGAGAAGACGCAAGACCAAGCGATGACGCGCTGCCCGCTCCACCCGCTCTCTGCGGTCCCGC
TGGTCCTTAGCCCCATGCCTCTTCCAAACAGTGCAATAAATCTCCGGTAACCCGTCCACCTGGCCGAG
GCAGCCAGAGCAGCGGGATGAGCTGGCGGCCAGAGAACGCCAGCCAGCCACCCCGCTCACGTCC
TGCCGACCCCTCCGCGGTCCCGACCTCTGCCCGCCAGCCCTCTGGCGTCTCTAGAACTGCTGCTG
ATGGAGACCAGGACAGAGCCCGCCCGCCCGCCAGCCCTCTGCCCTCTGCCAGGCCCTCC
CTCCCGGTGGTGGCGGCACCTCTGCTCCTCACGTGGGTGCGGGGCGGGCTTGTCAGGTTTGT
GAGACAGTCTTTTTTATTTGTATTGTTATTTTATTTTAAATTTAAACCTGATGACTTGACAG
CATCTTTCCACCGGGAAGAGCGGGCTGCTGCCTTCCCTCCTGCGGGGTGGGGTGGGACAGACCCA
GTGGGGGCCAGGGCCACACTCAGCGAGTGTGGGTGGACCCGCTCTGTGCTCCAGGGGTGGCCGGT
CTTCGGGGGCGCGGCTTCCACCACTCGCCGCCCGCCCGCTGCACTCCGTGACCTCTGAGAACCCAGCT
GGCCCTCTGTAGCCCCAGGTCTGCCGTCTAGAGGGAGCCACCGGCTCCGACACTGTCTTCCCGC
TAGAGCCCCGTCTCTCTGCGGGGCTCTCCGACCCCTTCCCTCACCTCCAGGCAGGGACAGAAT
AAATGTTTGTATGGATTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI



|                        |                                                                                                                                                                                                                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b> | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| <b>Components:</b>     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| <b>Note:</b>           | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                   |
| <b>RefSeq:</b>         | <u><a href="#">NM_020695.4</a></u>                                                                                                                                                                                                                                                         |
| <b>Summary:</b>        | Seems to have no detectable effect on transcription elongation in vitro.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                    |
| <b>Locus ID:</b>       | 57455                                                                                                                                                                                                                                                                                      |
| <b>MW:</b>             | 30.1                                                                                                                                                                                                                                                                                       |