

Product datasheet for **SC210165**

ZNF586 (NM_017652) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF586

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_017652

Insert Size: 850 bp

Insert Sequence: >SC210165 3'UTR clone of NM_017652
The sequence shown below is from the reference sequence of NM_017652. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTTCACTACTGGAATGAGGCCTTATAAGTGAAGCAAATTTTGAAAATCTCTTGCCAGGCTTTTGTCTC
CTTCAAGGCCAGAGAGTTACACCCAGATCAAGGTGTTATGAGTGTGACAAATGGGGAATATCTTTAGC
TAGAATGCTAGCTTCTTTACATAAAAGAGTGCTCCCACTGAAGAAGTGCCTTTTGAGTGCAATGAATGT
GAGAAAGCCTTCAGCCCTCTCTCATTGGTTACCACAATATTTACATGAGGAAAATACCACAGATGTGGA
GGTAATATTATTTTCTTCAATATAACACTGGAGGAAACCCCTTATGAGGATGCCATCTGCCTAAATT
GAATGTCATACATCAAATAGCTGCATACATTCCAGGTATGTGGAACTGCATAGTATTTTTCGACCTGC
CCAGGTCTCTTGCTAGACTTCTGTGACTTACAGTTCTCTTCTAAAAAGCATTTACCTTTACCACTTGG
CAGGTGCCACAGTACATCATTCGCCTAACTCATCATGATCCAGAAGTAACTTTGATTTTGTTCATTC
AACAGAGGGATTTTTAGTATCCTCAGCACTATTCTATTCTATTCTTTCTCTCCAATGAGTTAGGC
CATGAACCTGACTCATTTTTGGCACAGAGGACTCTTCTGATAACTTGAAAAATGGACTACCTTTTCA
GGGATCTTTTGGATCTCACAGAGTTCTTAAGATGGAATTTTCCAGCCTGCCAGTCATGAATCTCAGACA
GCCTGCCACCTATTGCCCTGATTTGTGTTATAATAAAGGCCTTGTGCTTAAGACACCTTCAGTCTTTG
CAGGAGGGCATGGAAGCTCAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_017652.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	54807
MW:	32.5