

Product datasheet for SC210167

NAPIL2 (NM_021963) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: NAPIL2

Synonyms: BPX

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_021963

Insert Size: 841 bp

Insert Sequence: >SC210167 3'UTR clone of NM_021963
The sequence shown below is from the reference sequence of NM_021963. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GAGGCATTAGTAGAAGACATTGATCGTAGAGCAGAGTATACATGGCCCTGAAATTAAGTGCCTAGAT
ATAGTTACTCAAGGTATAAGAAGCCTTGTTCTGTATTTTCTGTAGTGTAGTTAAACATATGTT
TCAAAAATATAAGAAAAGTTCAAAAATAATTAATTTGACCTTGAGTTTGTAGTAGAATGTTTCAA
GAAATGTACACTGTGGTAAATGATTTAAACACTAGTATAGTGTTGTAGCTTAATCCTTCTGAAGTC
TTTTTGTATGTAGCTATTAATCTGTGGCTATGAAATGATCAGAAATGCTAAGTGAGATCAATTTGT
TTGGAAAAAAATCTTGGGAAACAACCAAGGGTTTTTCGCTGTTGTTGTTTTCTTTTCTATTTTGT
TTACTTAGTCCTTAGCTAGTGGATTAATTTGTTGCTGCTTCATTTTGAATAACAATGCAGTA
GAATTTAAACTTTGGATGCTTAAGAGGCCTGCATATAGATAAGAATTCAGGCAAACTACATTTATTG
TTAATAACAGCTTGTTCATAGGCTCTTGTATTTATGTAAGTGTGATAAATAATGAACTTAGTTATAT
TGAGGTTATTGTTGTCGGTGAAGTGTAGTCACAGTATTTTCAAAAGTTTGCACATATTGTTCTGTGT
AATTGTGAAGCCATAATTACAGTGTTTAATTCCTTTTCTATTACATCATTGAAAGTGATCAC
TTTACCATTTTGAAAAGATATTTCTGTTCTTTCACTGCAAAATAAAAAGAATAAAATTTTCAGAGTGT
CTCATGTTTAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_021963.4</u>
Summary:	The protein encoded by this intronless gene is a member of the nucleosome assembly protein (NAP) family. The encoded protein represents a class of tissue-specific factors that interact with chromatin to regulate neuronal cell proliferation. [provided by RefSeq, Jan 2011]
Locus ID:	4674
MW:	32.6