

Product datasheet for **SC210185**

RAB5C (NM_004583) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RAB5C (NM_004583) Human 3' UTR Clone
Symbol:	RAB5C
Synonyms:	L1880; RAB5CL; RAB5L; RABL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004583
Insert Size:	834 bp
Insert Sequence:	>SC210185 3'UTR clone of NM_004583 The sequence shown below is from the reference sequence of NM_004583. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCAGCCGGAGCCAGTGTGCAGCAACTGAGCCCCCCTTGCCCTGCCCGCTGCCCCGCTCCTCCGCT
GAATGACCCGACTGGAATCCACTCTAACCAATCGCACTTAACGACTCGGGCCACCACTGGGGGGCAGG
GGGAGGGGTCCACCATGATTTCTCCATATAATTTTGATCATAGGCCGGAGTGAGTCATTCCACCTGCAC
CTTTCTGTACAAATACTAATTCAATTTTAAGTCTTAAGTCACTTTTTTAATATATATGATCTTCTGCTC
TTCCCACTTCTCCCTTTCTACTGCTCTCCATTTTCCCTTGCTGGGAGTAGCCACATGCTCTTGCCC
CCCAACCTTGTATATGGGGACAGTGGGGTCAGTGCAGCTACCCTTTCTTTCCCTCCTGCGGAACAGCG
GACCCAGCAAGAGCATCCACATCCTCACTTTGTTCCGGAGTGGTCTTTGGTTTGGGCGGTGGGGCAGACC
TTGGGAAGGGGCTTAGGAAGGGAGAGGCAGCTCTTCCCTCAGCTGGCTCTCATCAGGCTGCAGCCCCCT
CCCCGCTCCCACCTCCCTGCTGGGAAACCACAGCATTATCACAGCATTATTGTGACAGCCACGAACCCA
TTGCCACAACCCCTCCACCTCGGTACCCCCAACCTCTGGCTCTGAGCCCTGTTCTGACCAAATCATG
ATGATGAGTATTTGGGGTGGGTGGTAAGGGGGGAGTGGGAGGGGACGGAACCAACTTTTTCTGTAT
TTTGATTGTATGTTTTCTCAACATGTAACCAATCAGTATCTTGTCATATAGTCAGCCGATCGATCG
ACCTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004583.4</u>
Summary:	Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996 [PubMed 8646882]).[supplied by OMIM, Nov 2010]
Locus ID:	5878
MW:	30.5