

Product datasheet for **SC210193**

COMMD10 (NM_016144) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: COMMD10

Synonyms: PTD002

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_016144

Insert Size: 828 bp

Insert Sequence: >SC210193 3'UTR clone of NM_016144
The sequence shown below is from the reference sequence of NM_016144. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ATACAAGCACAGCTGGATTCCCTTACATGATGTTTTCGAAGACTGTTTTTTCATCACGCTCCTGCCAC
CTCATTATTTTGCATTGAAGATACATTGCCAGGTTGTGTTTTCTGAAGGATTCAGTGACTTGCTTTCTG
TAAATTATATGGCTTATCACTTCTTAGACAAATAACAACCAATAGAGATCATTGTTAAGAATACTGAGG
TTCTAATATACTTTCTTTAGTTCTGTGAGCCAACAGTAATTATTAAGAACACTTCCCTTTAAAGGAAA
CAAAAGTGAATACCATATTGTTTTACTGTCATAGTGTGCTTTCTTGCTGCTGCTGTAGTTTTTAC
TTGCTGGATGATACCATAATGTATCAAGGAGCGTCCATGGATACAAGATAAGATGTGTACCTTAGTAGA
ATACAGAGCTTTGGTAATTACATGAATAAAATTAAGAAAATAGCCATATACAATCAAATACACTATGGC
ATTTTTATTTGAATATGATGAGTATATTTTGCTTCGGAAATAATATAGGAAGGAAATGTAATAATAGTGA
GTAGTATGGTATCAGTTAATCCAGTCTGAGCTTCTCTGTCAACTTCAGTTTCTCTCAGTTTAATGA
TTTAATAATAGTCCAGGTTTTTGTGTGTTTTTCTTATACTGCAAATTAATAATGATTCACCTTTATAGT
TTGGGAGACAGAATCAGGTCTTGAATAAAATAATTGTAATGAGTGCTAAATGGGACCATTATTCGAAT
CAGATACCTTTTATATTCTCTTCCATAAATACGTTGATTTCTGTCAATAAAATTTTGTGTCTTAGGA
ACGCGTAAGCGGCCGCGGCATCTAGATTTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_016144.4
Summary:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966).[UniProtKB/Swiss-Prot Function]
Locus ID:	51397
MW:	32.3